



Helping all people
live healthy lives

Product Catalogue 2014

Eastern Europe, Middle East and Africa

BD Medical

Medical Surgical Systems



Contact Details

BD in Ghana

East and West Wings
4th Floor, Secant Tower,
7th Lane, Kuku-Hill, Osu – Accra
Greater Accra Region, Ghana
Phone +233 302 769692
Fax +233 302 769601

BD in Kenya

Africa Re-Center
Hospital Road, Upper Hill
PO Box 76613
Nairobi 00508, Kenya
Phone +254 2 02 738 339
Fax +254 2 02 738 342

BD in Russia

12, Building 1, Block C,
Dvintsev Street.,
127018 Moscow, Russia
Phone +7 495 775 8582
Fax +7 485 775 8583

BD in Turkey

Rüzgarlibahçe Mah. S. Sinan Eroglu Cad.
Akel Is Merkezi A Blok, Kat -3 No: 6
34805 Kavacik Beykoz - Istanbul, Turkey
Phone +90 216 680 10 02
Fax: +90 216 680 16 55

BD in South Africa

The Woodlands Office Park
Building 31, 2nd Floor
Woodlands Drive, Woodmead
PO Box 3569, Rivonia 2128
2199 Johannesburg, South Africa
Phone +27 11 603 2620
Fax +27 11 603 2640

BD in United Arab Emirates

Arif & Bint oak Building
1st floor
Karama, Zabeel Road,
P.O. Box 52279
Dubai, United Arab Emirates
Phone +971 43 379 525
Fax +971 43 379 551

BD in Zambia

Elunda Two, Second Floor
Addis Ababa Roundabout
Corner Chikwa Rd./Los Angeles Blvd.
Rhodes Park
Lusaka, Zambia
Phone +260 211 372 784.

Table of Contents

Injection Hypodermic Needles 5 Universal Fill Needles 7 Safety-Engineered Needles and Syringes 8 Needles and Syringes 12 Oral Dispensing Syringes 18 Syringes with Passive Reuse Prevention Technology 20 Autodisable Immunization Syringes 22	1
Infusion Therapy Extension Tubes 24 Infusion Sets 25 Closed Systems - Needleless IV Access 26 Stopcock for IV Therapy and Pressure Monitoring 28 Infusion Accessories 29 One-way Valve 30 Skin Prep 31	2
Peripheral Vascular Access Systems Closed IV Catheter Systems 32 Indwelling Cannulae 35 Indwelling Cannulae for Neonatal and Pediatric Use 43 IV Dressing 43 Obturator 44-45 BD Vialon™ Biomaterial 46 BD Cannulae / Catheter Conversions 47	3
Flushing of Vascular Access Devices BD PosiFlush™ Prefilled Flushing Syringes 48-49	4
Arterial Catheters Arterial Catheters 50-51 Arterial Cannula 50	5
Anaesthesia Systems Spinal Anaesthesia 52 Epidural Products 55 Combined Spinal / Epidural Anaesthesia Sets 57	6
Surgical Accessories Biopsy Needles 58 Heimlich Chest Drain 59 BD E-Z Scrub™ 59	7
Sharps Disposal Systems BD™ Sharps Disposal Systems 60-61	8
Closed System Drug Transfer Device BD PhaSeal™ System 62 Protector 63 Injector 64 Connector 65 Administration Products - IV Bag Access 66-67 Application Products 68	9
Social Responsibility Trusted Partners 72	10
Annex General Conditions of Sales 72-73 Quality Management with BD 74	11
Index Numerical Index 75 Alphabetical Index 79	12
Notes Notes 82	13



Injection Systems

Hypodermic Needles

BD Microlance™ 3

Easy penetration, smooth gliding force during insertion/withdraw and reduced cannula friction while sliding through tissue is a result of the following:

- Three bevel precision grinding of needle bevel
- Targeted needle surface polishing process
- Proprietary needle lubrication process

Thin walls:

Permitting the use of thinner needles with larger lumen thus increasing flow rates during injection and collection.

BD Microlance™ can be connected to a Luer adapter or Luer Lock connector. Adapter and color coding in compliance with ISO standards.

Recommended uses: intradermal, subcutaneous, intramuscular, intravenous injections and aspiration.

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Wall	Bevel	Color Code	Shelf Box / Shipping Case	Reference Number
BD Microlance™ Needle								
16 G	1,6	40	1 1/2"	regular	regular	white	100 / 5000	300637
18 G	1,2	40	1 1/2"	thin	short	pink	100 / 5000	304622
18 G	1,2	50	2"	regular	short	pink	100 / 4000	301900
19 G	1,1	25	1"	thin	regular	ivory	100 / 5000	301700
19 G	1,1	40	1 1/2"	thin	regular	ivory	100 / 5000	301500
19 G	1,1	50	2"	thin	regular	ivory	100 / 4000	301750
20 G	0,9	25	1"	thin	regular	yellow	100 / 5000	304827
20 G	0,9	40	1 1/2"	thin	regular	yellow	100 / 5000	301300
21 G	0,8	16	5/8"	thin	regular	green	100 / 5000	304434
21 G	0,8	25	1"	thin	regular	green	100 / 5000	301156
21 G	0,8	40	1 1/2"	thin	regular	green	100 / 5000	304432
21 G	0,8	50	2"	regular	regular	green	100 / 4000	301155
22 G	0,7	25	1"	thin	regular	black	100 / 5000	304727
22 G	0,7	30	1 1/4"	thin	regular	black	100 / 5000	300900
22 G	0,7	40	1 1/2"	thin	regular	black	100 / 5000	301000
22 G	0,7	50	2"	regular	regular	black	100 / 4000	300094
23 G	0,6	25	1"	thin	regular	blue	100 / 5000	300800
23 G	0,6	30	1 1/4"	thin	regular	blue	100 / 5000	300700
24 G	0,55	25	1"	regular	regular	violet	100 / 5000	304100
25 G	0,5	16	5/8"	regular	regular	orange	100 / 5000	300600
25 G	0,5	25	1"	regular	regular	orange	100 / 5000	300400
26 G	0,45	13	1/2"	regular	regular	brown	100 / 5000	303800
26 G	0,45	10	3/8"	regular	short	brown	100 / 5000	300300
26 G	0,45	16	5/8"	regular	regular	brown	100 / 5000	304300
27 G	0,4	13	1/2"	regular	regular	grey	100 / 5000	300635
27 G	0,4	19	3/4"	regular	regular	grey	100 / 5000	302200
30 G	0,3	13	1/2"	regular	regular	yellow	100 / 5000	304000

1

2

3

4

5

6

7

8

9

10

11

12

13

Injection Systems


Percutaneous Injuries

Would you believe over 1 million needlestick injuries (NSIs) occur every year¹?

"Injuries caused by needles and other sharp instruments are one of the most common and serious risks to healthcare workers and represent a high cost for health systems and society in general. According to healthcare experts, such injuries occur more than one million times per year, partly causing very serious consequences and sometimes leading to serious diseases."¹ Nurses and doctors are most at risk of needlestick injuries².

Facts:

- Injection and infusion devices cause more than 1/3 of sharps injuries²
- Any sharps injury could be a potential risk no matter the source³
- Not all hazardous exposures to healthcare workers (HCW) are percutaneous injuries²
- Some exposures are to the eye or other mucous membranes²
- A number of NSIs occur after use, during the disposal process and anyone in the hospital environment could be at risk (e.g. cleaning staff)²



A comprehensive safety approach shall have three tenants:



BD Injection Expertise. Our Experience in your Hands

¹EU Council Directive 2010/32/EU, June 11, 2010.

² <http://healthsystem.virginia.edu/internet/safetycenter/>

³ Wittmann, A, Köver, J, Hofmann, F, Kralj, N (2009): Übertragene Blutvolumina nach Nadelstichverletzungen an s.c. Kanülen. Dok.Dt.Ges.Arbmed.Umweltmed. 49. Jahrestagung: 382-84

Injection Systems

Universal Fill Needles

BD™ Blunt Fill Needle

- Fill Needle with blunt tip
- 45° bevel avoids coring and designed to minimize the risk of needle stick injury
- Red hub and shield for easy identification as fill needle
- Also available with 5 micron filter for effective filtration of glass, metal, rubber and other contaminant particles
- Sterilization method for
 - REF 305211 is Ethylene Oxide
 - REF 305180/181 is Radiation
- Sterile, single use
- CE marked



Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Bevel	Color Code	Shelf Box / Shipping Case	Reference Number
BD™ Blunt Fill Needle							
18 G	1,2	40	1 1/2"	blunt	red	100 / 1000	305180
18 G	1,2	25	1"	blunt	red	100 / 1000	305181
BD™ Blunt Fill Needle with 5 Micron Filter							
18 G	1,2	40	1 1/2"	blunt	lilac	100 / 1000	305211

1

2

3

4

5

6

7

8

9

10

11

12

13

Injection Systems

Safety-Engineered Needles and Syringes

BD Eclipse™ – The Safety Needle with SmartSlip™ Technology

The SmartSlip™ Technology is especially designed for safe assembling of BD Eclipse™ to Luer syringes.

- Bevel oriented to the safety cover to allow low angle injections
- Activates with one finger stroke - no hard surface required
- Minimal change in technique
- No excess dead space - minimal medication waste

Safety

- By assembling BD Eclipse™ to a Luer slip syringe, an audible click indicates a safely assembled needle
- Safety activation immediately follows injection

- Wide, textured finger pad for thumb or forefinger
- Single handed activation of the safety shield, audible click indicates that the cover has safely engaged the needle

Patient Care

- Activation occurs after needle is withdrawn
- BD needle technology is designed for patient comfort
- Sterilization method is Ethylene Oxide
- Sterile, single use
- Does not contain natural rubber latex
- CE marked



Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Eclipse™ Safety Needle						
30 G	0,3	13	1/2 "	yellow	100 / 1200	305771
27 G	0,4	13	1/2 "	grey	100 / 1200	305770
27 G	0,4	19	3/4 "	grey	100 / 1200	305889
27 G	0,4	40	1 1/2 "	grey	100 / 1200	302436
25 G	0,5	16	5/8 "	orange	100 / 1200	305760
25 G	0,5	25	1 "	orange	100 / 1200	305891
23 G	0,6	25	1 "	blue	100 / 1200	305892
23 G	0,6	30	1 1/4 "	blue	100 / 1200	305886
22 G	0,7	30	1 1/4 "	black	100 / 1200	305887
21 G	0,8	25	1 "	green	100 / 1200	305894
21 G	0,8	40	1 1/2 "	green	100 / 1200	305895
20 G	0,9	25	1 "	yellow	100 / 1200	305899
20 G	0,9	40	1 1/2 "	yellow	100 / 1200	305888
18 G	1,2	40	1 1/2 "	pink	100 / 1200	302437

Gauge Size	Capacity	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box	Reference Number
BD Eclipse™ Safety Needle with BD Luer-Lok™ Syringe							
25 G	3 ml	0,5	16	5/8 "	orange	50 / 300	305781
23 G	3 ml	0,6	25	1 "	blue	50 / 300	305782
22 G	3 ml	0,7	40	1 1/2 "	black	50 / 300	305783
21 G	3 ml	0,8	40	1 1/2 "	green	50 / 300	305784
22 G	5 ml	0,7	40	1 1/2 "	black	50 / 300	305785
22 G	10 ml	0,7	40	1 1/2 "	black	50 / 300	305786

Injection Systems

Safety-Engineered Needles and Syringes

BD SafetyGlide™ – Safety Needle with Activation-Assist™ Technology

BD SafetyGlide™ Safety Needles with patented Activation-Assist™ Technology help to protect from needle stick injuries through its integrated safety mechanism.

- Low profile safety mechanism remains discretely out of the way before and during injection
- Patented BD Activation Assist™ technology ensures quick coverage of the needle by a single finger stroke immediately after injection
- Minimal change in technique
- “Bevel up” design eliminates safety-arm interference to facilitate low angle injections
- Sterilization method is Radiation
- Sterile, single use
- Does not contain natural rubber latex
- CE marked



Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD SafetyGlide™ Safety Needle						
27 G	0,4	16	5/8"	grey	50 / 500	305921
25 G	0,5	16	5/8"	orange	50 / 500	305901
25 G	0,5	25	1"	orange	50 / 500	305916
23 G	0,6	25	1"	blue	50 / 500	305902
22 G	0,7	40	1 1/2"	black	50 / 500	305900
21 G	0,8	25	1"	green	50 / 500	305915
21 G	0,8	40	1 1/2"	green	50 / 500	305917
18 G	1,2	40	1 1/2"	pink	50 / 500	305918

BD SafetyGlide™ - with BD Luer-Lok™ Syringe packed aside

Safety

- Safety activation immediately follows injection
- Needle is safely protected after shielding
- Single-handed activation

Patient Care

- Activation occurs after needle is withdrawn
- Sterilization method is Radiation
- Single, sterile use
- Does not contain natural rubber latex
- CE marked



Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD SafetyGlide™ with BD Luer-Lok™ Syringe							
3 ml	25 G	0,5	16	5/8"	orange	50 / 400	305904
3 ml	23 G	0,6	25	1"	blue	50 / 400	305905
5 ml	22 G	0,7	40	1 1/2"	black	50 / 400	305907

Injection Systems

Safety-Engineered Needles and Syringes

BD SafetyGlide™ - Syringe for TB and Allergy

BD SafetyGlide™ syringes complement BD's safety injection family with the Tiny Needle Technology (TNT) to the most versatile safe injection technology.

The BD SafetyGlide™ Syringe for TB and Allergy is specially designed and scaled to deliver smaller doses of medication. It has an integrated needle and helps to meet the unique needs of the healthcare workers who administer these types of injections.

- Minimal change in technique
- Finger stays behind the needle at all times
- Secure and visible lock
- Rotating safety arm eliminates interference during drawing up and facilitates low angle injections
- Integrated needle means minimal waste space
- 26G and 27G needles perform a variety of subcutaneous and intradermal applications
- Sterilization method is Radiation
- Sterile, single use
- Does not contain natural rubber latex
- CE marked



Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD SafetyGlide™ with integrated needle							
1 ml	27 G	0,4	13	1/2"	grey	100 / 400	305945
1 ml	26 G	0,45	10	3/8"	brown	100 / 400	305946

Injection Systems

Safety-Engineered Needles and Syringes

BD Integra™ - with Detachable, Retracting Needle

BD Integra™ syringe with retracting needle combines innovative technology with BD quality and comfort.

- Spring based retracting syringe with detachable needle
- Unique non-luer thread design - BD Tru-Lok™ of the BD Integra™ detachable needle provides:
 - secure connection and
 - is designed to have significantly lower waste space volume as compared to a traditional luer lock syringe and needle. Less waste space can reduce drug wastage
- Activation of retracting needle before or after needle is withdrawn from the patient - clinician takes the decision
- Complete dosing ensured prior to activation
- BD Integra™ blunt needle allows change of needle for aspiration
- No change in technique

Safety

- Safety activation immediately follows injection
- Needle is completely retracted into the syringe barrel

Patient Care

- Syringe and needle cannot be reused
- BD PrecisionGlide™* needle technology is designed for patient comfort
- Sterilization method for
 - REF 305833 is Ethylene Oxide
 - others it is Radiation
- Sterile, single use
- Does not contain natural rubber latex
- CE marked

** BD PrecisionGlide™ is BD's US brand for hypodermic needles corresponding to the European brand BD Microlance™*



Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Integra™ syringe only							
3 ml	–	–	–	–	–	100 / 400	305283
BD Integra™ with retracting needle							
–	25 G	0,5	16	5/8"	orange	100 / 1000	305310
–	25 G	0,5	25	1"	orange	100 / 1000	305311
–	23 G	0,6	25	1"	blue	100 / 1000	305312
–	22 G	0,7	40	1 1/2"	black	100 / 1000	305313
BD Integra™ Blunt Fill Needle (non-retracting needle)							
–	18 G	1,2	40	1 1/2"	red	100 / 1000	305275
BD Integra™ Blunt Fill Needle with 5 micron filter (non-retracting needle)							
–	18 G	1,2	40	1 1/2"	red	50 / 500	305833
BD Integra™ Syringe with mounted, retracting BD PrecisionGlide™ needle							
3 ml	25 G	0,5	16	5/8"	orange	100 / 400	305269
3 ml	25 G	0,5	25	1"	orange	100 / 400	305270
3 ml	23 G	0,6	25	1"	blue	100 / 400	305271
3 ml	22 G	0,7	40	1 1/2"	black	100 / 400	305272
3 ml	21 G	0,8	25	1"	green	100 / 400	305273
3 ml	21 G	0,8	40	1 1/2"	green	100 / 400	305274

Injection Systems

Needles and Syringes

BD Emerald™ Syringe For General Purpose Use

Combines High-Quality Performance and Reduced Environmental Impact



A new general purpose syringe inspired by BD's ongoing commitment to global health

Achieving Global Performance Standards

Every second of every day, approximately 220 BD syringes are used worldwide. Based on this global presence, as well as over 100 years of hypodermic innovation, BD has created design and process specifications for the BD Emerald™ Syringe that help ensure it will consistently meet performance standards such as ISO 7886-1 and ISO 594.

BD Emerald™ Syringes must meet over forty specific criteria before they are qualified for

use. BD's passion for ongoing improvement continues to raise the bar of product reliability.

Reducing Environmental Impact

Because human health and a healthy environment are inseparable, BD is dedicated to the design, manufacture and distribution of safe, effective medical technology products which help to minimize environmental impact.

The BD Emerald™ product line has been produced in a more

sustainable manner, with up to 30% reduction in material.

If half of the world's syringes could be converted to BD Emerald™ syringes, 15,000,000kg of waste could be saved per year. This is the equivalent combined weight of approximately 200,000 people – enough to fill two Olympic stadiums!

BD Emerald™ Syringe.

A Step Toward a Greener Future.

Injection Systems

Needles and Syringes

BD Emerald™ – 2, 5, 10 ml

BD Emerald™ syringe for general purpose use. Specialized design features benefit clinicians, patients and the environment.

- Green stopper and bold scale markings to aid in setting an accurate dose
- Clear barrel to allow for visualization of fluid
- Robust retaining ring to help prevent accidental plunger rod pullout
- Flange size and shape to provide stability and comfort during use
- Textured thumb press to reduce slippage during administration
- Meets international standards: ISO 7886-1 and ISO 594 compliant
- Minimizes waste disposal – has up to 30% less material than other syringes*

- Sterilization method is Ethylene Oxide
- Sterile, single use
- Does not contain natural rubber latex
- CE marked

Each BD Emerald™ Syringe:

- Addresses clinical needs globally
- Meets quality standards throughout the world
- Reduces environmental impact
- PLUS - Offers you BD's full array of World Class Sharpness needles



* Data on file. Based on comparison with leading syringe brands worldwide

Capacity	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Emerald™ syringe				
2 ml Luer Slip	concentric	0,1 ml	100 / 3000	307727
5 ml Luer Slip	concentric	0,2 ml	100 / 2000	307731
10 ml Luer Slip	concentric	0,2 ml	100 / 1200	307736

Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Scale Graduation	Color Code	Shelf Box / Shipping Case	Reference Number
BD Emerald™ syringe with pre-attached needle								
2 ml Luer Slip	22 G	0,7	30	1 1/4"	0,1 ml	black	100 / 2000	307728
2 ml Luer Slip	23 G	0,6	25	1"	0,1 ml	blue	100 / 2000	307740
2 ml Luer Slip	23 G	0,6	30	1 1/4"	0,1 ml	blue	100 / 2000	307741
5 ml Luer Slip	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 1500	307732
5 ml Luer Slip	22 G	0,7	30	1 1/4"	0,2 ml	black	100 / 1500	307733
5 ml Luer Slip	23 G	0,6	30	1 1/4"	0,2 ml	blue	100 / 1500	307735
10 ml Luer Slip	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 900	307737
10 ml Luer Slip	22 G	0,7	30	1 1/4"	0,2 ml	black	100 / 900	307738

Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Scale Graduation	Color Code	Shelf Box / Shipping Case	Reference Number
BD Emerald™ syringe with detached needle								
2 ml Luer Slip	23 G	0,6	25	1"	0,1 ml	blue	100 / 3000	307729
2 ml Luer Slip	23 G	0,6	30	1 1/4"	0,1 ml	blue	100 / 3000	307730
5 ml Luer Slip	22 G	0,7	40	1 1/2"	0,2 ml	black	100 / 2000	307734
5 ml Luer Slip	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 2000	307742
5 ml Luer Slip	22 G	0,7	30	1 1/4"	0,2 ml	black	100 / 2000	307743
10 ml Luer Slip	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 1200	307739

Injection Systems

Needles and Syringes

BD Discardit™ II - 2-Piece Syringe

- Transparent barrel for good visualization of syringe contents
- Smooth plunger movement - reduced sliding force
- BD Discardit™ syringes are tested during manufacture in compliance with ISO 7886-1 to make them leak-tight
- Accurate graduation suited to user's needs
- With a retaining ring to prevent accidental plunger rod withdrawal
- Design of finger grip flanges ensures stability and comfort during injection
- Ergonomically designed plunger for single-handed operation
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Discardit™ 2-piece syringe				
2 ml Luer Slip	concentric	0,1 ml	100 / 3000	300928
5 ml Luer Slip	eccentric	0,2 ml	100 / 1800	309050
10 ml Luer Slip	eccentric	0,5 ml	100 / 1200	309110
20 ml Luer Slip	eccentric	1 ml	80 / 960	300296

Capacity	Gauge Size	Outer diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Discardit™ 2-piece syringe with BD Microlance™ needle packed aside							
2 ml Luer Slip	21G	0,8	40	1 1/2"	green	100 / 3000	300318
2 ml Luer Slip	22G	0,7	30	1 1/4"	black	100 / 3000	300320
2 ml Luer Slip	22G	0,7	40	1 1/2"	black	100 / 3000	300321
2 ml Luer Slip	23G	0,6	25	1"	blue	100 / 3000	300322
2 ml Luer Slip	23G	0,6	30	1 1/4"	blue	100 / 3000	300323
5 ml Luer Slip	22G	0,7	40	1 1/2"	black	100 / 1800	300506
5 ml Luer Slip	21G	0,8	40	1 1/2"	green	100 / 1800	300512
5 ml Luer Slip	22G	0,7	30	1 1/4"	black	100 / 1800	300513
5 ml Luer Slip	23G	0,6	30	1 1/4"	blue	100 / 1800	300516
10 ml Luer Slip	22G	0,7	30	1 1/4"	black	100 / 1200	300626
10 ml Luer Slip	21G	0,8	40	1 1/2"	green	100 / 1200	301001
10 ml Luer Slip	22G	0,7	40	1 1/2"	black	100 / 1200	301504
20 ml Luer Slip	21G	0,8	40	1 1/2"	green	80 / 960	300330

Injection Systems

Needles and Syringes

BD Plastipak™ - 2, 5, 10 ml

- Latex free stopper: eliminates completely the risk of allergic reactions to latex
- Clear barrel for visualization of syringe contents
- Bold scale printing to aid in setting for a more accurate dosage of medication
- Special silicone lubrication ensures that plunger moves smoothly and evenly
- Retaining ring: prevents accidental plunger rod withdrawal
- Sterilization method is Ethylene Oxide
- Sterile, single use
- Does not contain natural rubber latex
- CE marked
- Not autoclavable



Capacity	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Plastipak™ luer-slip syringe				
2 / 2,5 ml	concentric	0,1 ml	100 / 800	300185
5 ml	concentric	0,2 ml	100 / 400	302187
10 ml	eccentric	0,5 ml	100 / 400	302188

BD Syringe Luer-Lok™ Tip - 3, 5, 10 ml

- Clear barrel for visualization of syringe contents
- Bold scale printing to aid in setting for a more accurate dosage of medication
- Special silicone lubrication ensures that plunger moves smoothly and evenly
- Retaining ring: prevents accidental plunger rod withdrawal
- Leak-tight (locking facility on BD Luer-Lok™ syringes)
- BD Luer-Lok™ - needle hub or extension line can be safely locked on the syringe
- Sterilization method is Radiation
- Sterile, single use
- Does not contain natural rubber latex
- CE marked
- Not autoclavable

Capacity	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Syringe Luer-Lok™ tip				
3 ml	concentric	0,1 ml	200 / 800	309658
5 ml	concentric	0,2 ml	125 / 500	309649
10 ml	concentric	0,2 ml	100 / 400	300912

1

2

3

4

5

6

7

8

9

10

11

12

13

Injection Systems

Needles and Syringes

Low Dose Syringes

1 ml syringes make it easier to accurately measure a small dose than with a large capacity syringe. Available with or without needle.

BD Plastipak™ 1 ml Sub-Q and 1 ml Test Syringes with integrated needle

- Scale graduation 0,01 ml
- Without dead space
- Minimizes the formation of air bubbles
- Prevents any risk of detachment of needle
- Ensures administration of full dose
- BD Plastipak™ "Sub-Q" syringe with a short bevel needle for intradermal injections
- Packed in polybags with ten syringes each
- Sterilization method is Radiation

BD Plastipak™ 1 ml Syringes with mounted needle

- Compact "syringe plus needle" unit provides all the benefits of a BD Plastipak™ syringe and a BD Microlance™ needle at the same time
- Pre-assembled syringe plus needle in a single pack simplifies injections by reducing the number of steps and thereby reduces costs
- Sterilization method is Ethylene Oxide
- Sterile, single use
- Does not contain natural rubber latex
- CE marked



Description	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Plastipak™ 3-piece syringe				
1 ml Tuberculin / Allergy	concentric	0,01 ml	100 / 800	300013
1 ml Insulin Syringe U-40	concentric	U-40	100 / 800	300026
1 ml Insulin Syringe U-100	concentric	U-100	100 / 800	301355

Description	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Luer-Lok™ syringe				
1 ml BD Luer-Lok™	concentric	0,01 ml	100 / 800	309628

Description	Gauge Size	Outer Diameter (mm)	Length (mm)	Length (inch)	Bevel	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Plastipak™ syringe with mounted BD Microlance™ needle								
1 ml	25 G	0,5	16	5/8"	regular	0,01 ml	100 / 800	300014
1 ml	26 G	0,45	10	3/8"	intradermal	0,01 ml	100 / 800	300015
1 ml Insulin Syringe U-40	30 G	0,3	13	1/2"	regular	U-40	100 / 800	300334

Description	Gauge Size	Outer Diameter (mm)	Length (mm)	Length (inch)	Bevel	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Plastipak™ '1 ml Sub-Q' Syringe and BD Plastipak™ '1 ml Test' syringe with integrated needle, Packed in polybags with ten syringes each*								
1 ml Sub-Q, brown	26 G	0,45	13	1/2"	intradermal	0,01 ml	100 / 500	305501
1 ml Test, grey	27 G	0,4	10	3/8"	regular	0,01 ml	100 / 500	305502

* The interior of each syringe remains sterile until the two caps on the syringe - on the needle and the plunger rod end - are either removed or damaged.

Injection Systems

Needles and Syringes

BD Large Size Syringes

- Leak-tight (ensuring a secure connection on BD Luer-Lok™ syringes)
- Gradual and controlled movement of plunger
- Needle hub or extension line can be safely locked on the syringe
- Designed for accuracy and safety of use: transparent body
- Eccentric tip allows for low angle administration and aspiration: reduces dead space when using large capacity syringe
- Catheter tip: permits attachment of a probe or other tubing
- 100 ml syringe: catheter tip with potential to change into luer tip via an adapter
- 50 ml amber syringe for guaranteed protection of photosensitive drugs: stops 90% of light rays at 400 nm UV wave length
- Luer lock syringes are particularly suitable for automatic infusion equipment: syringe-drivers, pumps, etc. They have been approved for use by leading brands of syringe-drivers
- Sterilization method for
 - REF 305501 is Ethylene Oxide
 - REF 305502 is Radiation
- Single, sterile use
- Does not contain natural rubber latex
- CE marked
- Not autoclavable



Description	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
-------------	-----	------------------	---------------------------	------------------

BD Plastipak™ 20, 30 and 100 ml syringe

20 ml Luer Slip	eccentric	1 ml	120 / 480	300613
20 ml Luer-Lok™	concentric	1 ml	120 / 480	300629
30 ml Luer Slip	eccentric	1 ml	60 / 240	301231
30 ml Luer-Lok™	concentric	1 ml	60 / 240	301229
100 ml Catheter Tip with luer slip adapter	concentric	2 ml	25 / 50	300605

Description	Gauge Size	Length (mm)	Length (inch)	Scale Graduation	Shelf Box / Shipping Case	Reference Number
-------------	------------	-------------	---------------	------------------	---------------------------	------------------

BD Plastipak™ 20 ml syringe with BD Microlance™ needle packed aside

20 ml Luer Slip	21 G	40	1 1/2"	1 ml	60 / 240	301190
-----------------	------	----	--------	------	----------	--------

Description	Tip	Scale Graduation)	Shelf Box / Shipping Case	Reference Number
-------------	-----	-------------------	---------------------------	------------------

BD Plastipak™ and BD Perfusion 50 ml syringe

50 ml BD Perfusion Luer-Lok™	concentric	1 ml	50 / 100	300137
50 ml BD Plastipak™ Luer-Lok™	concentric	1 ml	60 / 240	300865
50 ml BD Plastipak™ Luer-Slip	eccentric	1 ml	60 / 240	300866
50 ml BD Plastipak™	catheter tip	1 ml	60 / 240	300867

BD Perfusion 50 ml syringe, opaque

50 ml BD Perfusion Luer-Lok™, amber	concentric	1 ml	50 / 100	300139
50 ml BD Plastipak™ Luer-Lok™, amber	concentric	1 ml	60 / 240	300869

Description	Gauge Size	Length (mm)	Length (inch)	Scale Graduation	Shelf Box / Shipping Case	Reference Number
-------------	------------	-------------	---------------	------------------	---------------------------	------------------

BD Perfusion 50 ml Syringe with mounted needle

50 ml BD Perfusion Luer-Lok™	14 G	30	1 1/4"	1 ml	50 / 100	300136
------------------------------	------	----	--------	------	----------	--------

BD Perfusion 50 ml Syringe with mounted needle, opaque

50 ml BD Perfusion Luer-Lok™, amber	14 G	30	1 1/4"	1 ml	50 / 100	300138
-------------------------------------	------	----	--------	------	----------	--------

Injection Systems

Oral Dispensing Syringes

BD Oral Dispensing Syringes - Clear and Amber

BD Oral syringes with an oversized (non-luer) tip to prevent connection to IV access lines and mounting of hypodermic needles, thus reducing the risk of parenteral administration of oral medications.

- Amber colored syringe protects light sensitive drugs
- BD oral syringes are supplied with leak-tight cap to prevent loss of the medication contained in the syringe after preparation

Recommended uses

- Specially designed for oral administration of drugs or fluids in drinkable form
- BD oral syringes provide the essential safety needed for the central preparation of unit doses
- Paediatric and geriatric departments
- Non-sterile, single use
- Packed in polybags with 100 syringes each
- CE marked



Capacity	Scale Graduation	Polybag / Shipping Case	Reference Number
BD oral dispensing syringe, clear with tip cap			
1 ml	0,02 ml	100 / 500	305217
3 ml	0,1 ml	100 / 500	305220
5 ml	0,2 ml	100 / 500	305218
10 ml	0,2 ml	100 / 500	305219
BD oral dispensing syringe, amber with tip cap			
1 ml, amber	0,02 ml	100 / 500	305207
3 ml, amber	0,1 ml	100 / 500	305210
5 ml, amber	0,2 ml	100 / 500	305208
10 ml, amber	0,2 ml	100 / 500	305209

Injection Systems

Reuse Prevention

The Situation...

16 billion injections are administered worldwide every year. In the developing world, as many as 40 % of these are given with reused syringes - exposing those seeking medical care to possible infection with the viruses that cause HIV/AIDS, Hepatitis B and Hepatitis C.



Proportion of injections given with reused equipment¹



Some syringe reuse is intentional due to economics or lack of proper training and education. Other reuse is unwitting: in some regions of the world, used syringes are repackaged and resold - appearing to be brand-new to the unsuspecting healthcare

worker. This practice spreads bloodborne diseases such as HIV/AIDS, Hepatitis B and Hepatitis C. The impact of these unsafe injections in terms of lives lost and destroyed is enormous.

Burden of Disease Due to Reused Injection Syringes in 2000²

	infections due to reused equipment	percent of disease attributable to reused equipment	healthy live years lost
Hepatitis B	22 million	32 %	3.1 million
Hepatitis C	2 million	40 %	300,000
HIV / AIDS	260,000	5 %	5.7 million

¹ World Health Organization (WHO) and Centers for Disease Control (CDC), 2002. Waste image provided by WHO

² WHO Global Burden of Disease Study, 2000, as cited in "A Procurement Guide: Procuring Single-use Injection Equipment and Safety Boxes", April 24, 2003, WHO

Injection Systems

Syringe With Passive Reuse Prevention Technology for Curative Use

BD Emerald™ PRO Reuse Prevention Syringe



Designed with Both Clinicians and Patients In Mind

BD Emerald™ PRO was developed specifically for curative procedures, which account for the vast majority of syringe injections.

The design features are the result of input from clinicians around the world. Opinions came from key countries where the BD Emerald™ PRO Syringe will help to protect patients from cross contamination of infections.

The resulting design is similar to a conventional syringe, yet features a strong¹ lock that meets the International standard and automatically renders the syringe unusable after an injection has been completed.

The BD Emerald™ PRO Syringe is easy and intuitive to use. It is designed with a green stopper and bold scale markings to aid in setting an accurate dose. In a clinical trial study,

clinicians felt it required only a minimal change in technique and provided both good control with smooth plunger movement for comfort and handling.

That's why BD Emerald™ PRO Syringe can be the solution for safe injection practices that you have been searching for.

¹ per ISO 7886-4

Injection Systems

Syringe With Passive Reuse Prevention Technology for Curative Use

BD Emerald™ PRO – Passive Reuse Prevention Syringe

The BD Emerald™ PRO syringe was designed to be flexible in use, with a robust locking mechanism. The resulting design is similar to a conventional syringe, yet features a strong¹ lock that meets the international standard and automatically renders the syringe unusable after an injection has been completed.

Other features are:

- Green stopper and bold scale markings to aid in setting an accurate dose
- Visual indicator to identify it as a reuse prevention syringe
- Passive reuse prevention feature to automatically lock after procedure is completed
- Breakable plunger rod to deter reuse after syringe is locked for added security
- Robust retaining ring to help prevent accidental plunger rod pullout
- Contains up to 30 percent less material than other syringes*
- Reduces environmental impact through reduction of material use

- Sterilization method is Ethylene Oxide
- Single, sterile use
- Does not contain natural rubber latex
- CE marked

Each BD Emerald™ Syringe:

- Addresses clinical needs globally
- Meets quality standards throughout the world
- Reduces environmental impact
- PLUS - Offers you BD's full array of World Class Sharpness needles



¹ per ISO 7886-4

* Data on file. Based on comparison with leading syringe brands worldwide

Capacity	Tip	Scale Graduation	Shelf Box / Shipping Case	Reference Number
BD Emerald™ PRO syringe				
2 ml	concentric	0,1 ml	100 / 3000	307960
5 ml	concentric	0,2 ml	100 / 2000	307970
10 ml	concentric	0,2 ml	100 / 1200	307980

Capacity	Gauge Size	Outer Diameter (mm)	Length (mm)	Length (inch)	Scale Graduation	Color Code	Shelf Box / Shipping Case	Reference Number
BD Emerald™ PRO syringe with needle packed aside								
2 ml	21 G	0,8	25	1"	0,1 ml	green	100 / 3000	307961
2 ml	21 G	0,8	40	1 1/2"	0,1 ml	green	100 / 3000	307962
2 ml	23 G	0,6	25	1"	0,1 ml	blue	100 / 3000	307964
5 ml	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 2000	307971
10 ml	21 G	0,8	40	1 1/2"	0,2 ml	green	100 / 1200	307981

Injection Systems

Autodisable Immunization Syringes

BD SoloShot™ Mini – Designed with Children and the Earth in Mind

- Easy to use, no change to clinical technique is required*
- Designed with minimal dead space to reduce medication waste*
- Ideal for large-scale and routine immunization programs
- The same quality and dependability that you've come to expect from BD for both intramuscular and intradermal injections
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- WHO PQS approved

** all claims are supported by data on file at Becton Dickinson, unless stated otherwise*




Capacity	Gauge Size	Outer Diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf Box / Shipping Case	Reference Number
BD SoloShot™ Mini, autodisable syringe with integrated needle							
0,5 ml	23 G	0,6	25	1"	blue	100 / 3300	302245
0,5 ml	24 G	0,55	19	3/4"	purple	100 / 3300	302247
0,5 ml	25 G	0,5	25	1"	orange	100 / 3300	302248
0,5 ml	25 G	0,5	16	5/8"	orange	100 / 3300	302249
BD SoloShot™ Mini, autodisable syringe with integrated needle - for low dose vaccine administrations							
0,05 ml	27 G	0,4	10	3/8"	grey	100 / 3300	302250
0,1 ml	27 G	0,4	10	3/8"	grey	100 / 3300	302251

New Mini Syringe is smaller than similar injection devices

The new compact size - smaller than similar injection devices - is designed to reduce the syringe's impact on the environment

- reduces use of raw material such as plastic and metal
- offers 38 % more syringes per case than similar devices, thereby reducing costs associated with shipping and transport
- provides storage efficiencies
- reduces disposal impact by virtue of its smaller size



Injection Systems

Safety Box For Disposal of Hypodermic Needles

BD Hub Cutter

- Portable, disposable, point-of-use needle hub destruction device
- Designed to cut the syringe needle at the hub and enclose the needle securely in the hard plastic container
- Cuts and holds 400-600 needles, depending on needle sizes
- Ready for use
- Separates needles and syringes at the point of use
- Large opening and surrounding wall facilitate needle insertion
- Compact and lightweight design
- Puncture-resistant container
- Designated fill line
- Yellow color coding for biohazardous waste
- Container cap can be closed between uses
- Container cap can be permanently locked prior to disposal
- Disposal of entire container avoids contact with sharps
- Easy to dispose - can be incinerated or buried in a pit
- Compatible with BD SoloShot™ Mini, BD Emerald™ PRO and any conventional syringe up to 50 ml



Description	Fill Capacity	Shelf Box / Shipping Case	Reference Number
BD™ Hub Cutter			
BD™ Hub Cutter for hypodermic needles	in average 400 - 600 needles	32	305307
BD™ Hub Cutter for hypodermic needles (Russian label)	in average 400 - 600 needles	32	305309

1

2

3

4

5

6

7

8

9

10

11

12

13

Infusion Therapy

Extension Tubes

BD Connecta™ Extension Tube to use with syringe pumps

- To use with syringe pumps
- Tube material: Phthalate (DEHP)-free PVC, transparent
- With BD Luer-Lok™ tip male / female
- Pressure resistant up to 2,0 bar

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Material	Priming Volume (ml)	Length (cm)	Outer Diameter (mm)	Inner Diameter (mm)	Shelf Box / Shipping Case	Reference Number
BD Connecta™ extension tube						
PVC	0,3	30	2,45	1,15	100 / 100	397994
PVC	1,55	150	2,45	1,15	100 / 100	397608
PVC	2,1	200	2,45	1,15	50 / 50	397610
BD Connecta™ Coex extension tube						
PE/PVC	1,55	150	2,45	1,15	100 / 100	397611

BD Connecta™ - Extension Tube for Infusion

- Tube material: Phthalate (DEHP)-free PVC, transparent
- Pressure resistant up to 2,0 bar

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Priming Volume (ml)	Length (cm)	Outer Diameter (mm)	Inner Diameter (mm)	Shelf Box / Shipping Case	Reference Number
BD Connecta™ extension tube for infusion					
1,75	25	3,5	3,0	50 / 50	397251
3,50	50	3,5	3,0	50 / 50	397509
5,30	75	3,5	3,0	50 / 50	397997
7,00	100	3,5	3,0	50 / 50	397590
10,50	150	3,5	3,0	50 / 50	397400

Infusion Therapy

Infusion Sets

BD Infusion Sets - Overview

- BD has a range of quality infusion sets which offer a maximum of safety and user comfort
- BD R 87 infusion contributes to managing risks within IV therapy due to the lockable ventilation and precise roller clamp
- The 15 mm liquid filter in the drip chamber as well as the BD Luer-Lok™ connector provide additional patient safety
- BD R 87 RLS infusion sets with back-valve integrated in the BD Luer-Lok™ offer an effective protection for parallel infusions
- BD R 87 Flexinert™ P infusion sets - the PVC-free alternative
- BD R 87 V infusion sets with automatic ventilation offer additional safety when there is risk for moistening the filter resulting in a halted infusion
- This is prevented by an integrated back flow valve in front of the bacterial tight filter in the air channel
- In addition we can also offer customized solutions
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Tube Length (cm)	Material	G	P	BD Luer-Lok™ Tip	Automatic Ventilation	Back Flow Valve	Injection Membrane	Shelf Box / Shipping Case	Reference Number
BD Infusion Sets Overview										
R 87 G with precise roller clamp	180	PVC	✓	–	✓	✓	–	–	100 / 100	396361
R 87 G with precise roller clamp	180	PVC	✓	–	✓	✓ lockable	–	–	100 / 100	396392
R 87 G with precise roller clamp	180	PVC	✓	–	✓	✓	✓	–	100 / 100	396370
R 87 P	150	PVC	–	✓	✓	–	–	–	100 / 100	396351
R 87 P	150	PVC	–	✓	✓	–	✓	–	100 / 100	396359
R 87 P	180	PVC	–	✓	✓	–	–	–	100 / 100	396350
R 87 P	180	PVC	–	✓	✓	–	✓	–	100 / 100	396353
R 87 P	220	PVC	–	✓	✓	–	–	–	100 / 100	396390
R 87 P with integr. 3-way stopcock and 10 cm tube	180	PVC	–	✓	✓	–	–	–	80 / 80	396354
R 87 P with integr. 3-way stopcock and 7 cm tube	180	PVC	–	✓	✓	–	–	–	80 / 80	396357
R 87 P Flexinert™	180	PVC free	–	✓	✓	–	–	–	100 / 100	396352
R 87 V	150	PVC	–	✓	✓	✓	–	–	100 / 100	396368
R 87 V	180	PVC	–	✓	✓	✓	–	–	100 / 100	396363
R 87 V	180	PVC	–	✓	✓	✓	✓	–	100 / 100	396373
R 87 V	180	PVC	–	✓	✓	✓ lockable	–	–	100 / 100	396371
R 87 V	180	PVC	–	✓	✓	✓ lockable	✓	–	100 / 100	396377

G = Gravity , P = Pressure

Infusion Therapy

Closed Systems - Needleless IV Access

BD Q-Syte™ - Luer Access Split Septum

BD Q-Syte™ Luer Access Split Septum is a needleless IV access device which is uniquely split. Split-septum features such as simple internal design, ease of use and a straight, clear fluid path are critical to achieve catheter-related bloodstream infection (CRBSI) reductions.

BD extended the benefits of split septum to the convenience of luer access with BD Q-Syte™ Luer Access Split Septum. In addition to the benefits of split septum design, the BD Q-Syte™ device delivers optimal luer access performance.

Because of its straight and unobstructed fluid path, the BD Q-Syte™ device provides dramatically higher flow rates¹, a low priming volume and the flexibility to use ISO 594 compatible Luer-Slip or Luer Lock connections

Features of BD Q-Syte™ Luer Access Split Septum

- Smooth silicone surface
 - swabable, easily disinfected
 - no interstitial space
 - may help to achieve appropriate disinfection

- Split Septum
 - Luer activated closed system
 - helps prevent micro-organism ingress
- Lipid compatibility
 - suitable for parenteral nutrition
- Luer lock and luer slip compatibility
- Straight large fluid path
 - allows high flow rates (> 500 ml/min)
- Clear housing
 - complete internal visibility
 - helps to assure flush efficiency
- Compact light device
 - may help to increase patient comfort due to the smaller and lighter design
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- Reference Number 385102 contains DEHP

¹ data on file



Description	Number of BD Q-Syte™	Flow Rate HzO (l/hr)	Priming Volume (ml)	Shelf Box / Shipping Case	Reference Number
BD Q-Syte™ Luer Access Split Septum					
BD Q-Syte™ Stand-alone device	–	32	0,16	50 / 200	385100
BD Q-Syte™ Extension Sets					
BD Q-Syte™ 15 cm (6 in.) Standard Bore, fixed nut	1	27	1,14	25 / 200	385101
BD Q-Syte™ 15 cm (6 in.) Small Bore, fixed nut	1	3	0,34	50 / 200	385102*
Acacia Inc.** Specialty Extension Sets with BD Q-Syte™ Luer Access Split Septum					
15 cm (6 in.) Macro Bore Bi-Extension Set, spin nut	2	20	1,60	50 / 50	385161
15 cm (6 in.) Macro Bore Tri-Extension Set, spin nut	3	20	2,25	50 / 50	385162
15 cm (6 in.) Small Bore Bi-Extension Set, spin nut	2	2	0,45	50 / 50	385163
15 cm (6 in.) Small Bore Tri-Extension Set, spin nut	3	2	0,80	50 / 50	385164

* contains DEHP

** Acacia Inc. product featuring the BD Q-Syte™ Luer Access Split-Septum are exclusively distributed by BD Medical. Acacia Inc., Brea, CA. 92921 www.acaciainc.com

Infusion Therapy

Closed Systems - Needleless IV Access

BD Q-Syte™ - Luer Access Split Septum

Closed System

- Use on ISO standard IV medical devices
- Use with intermittent or continuous infusion
- Connection and disconnection of syringes and IV sets
- Blood withdrawal
- Use on arterial lines




Description	Number of BD Q-Syte™	Priming Volume (ml)	Shelf Box / Shipping Case	Reference Number
BD Connecta™ 3-way with BD Q-Syte™				
BD Connecta™ Stopcock with BD Q-Syte™	1	0,31	50 / 500	394501
BD Q-Syte™ Accessory				
BD Q-Syte™ Vial Access Adapter	–	–	25 / 100	385108

Based on peer-reviewed evidence, the CDC¹ has revised the Guidelines for the Prevention of Intravascular catheter-related bloodstream infections to include a recommendation on the preferred type of needleless connector used.


While many needleless connectors may claim to be a split septum, the split septum devices described in the evidence cited by the CDC were BD Q-Syte™ Luer Access Split Septum and BD Interlink®.

A split septum needleless access system has 64% - 70% lower CRBSI rates than mechanical valves.^{2,3}



BD Q-Syte™ Luer Access Split Septum

- no internal components or mechanisms
- straight fluid path
- smooth connection surface with clear external housing



Mechanical valve needleless connector

- indirect fluid path
- irregular gaps
- complex internal components
- internal fluid reservoirs

¹ CDC Guidelines May 2011

² Rupp ME, Sholtz LA, Jourdan DR et al. (2007): Outbreak of bloodstream infection temporally associated with the use of an intravascular needleless valve. *Clin Infect Dis.* 44(11):1408-14.

³ Salgado CD, Chinnes L, Paczesny TH, Cantey JR (2007): Increased rate of catheter-related bloodstream infection associated with use of a needleless mechanical valve device at a long-term acute care hospital. *Infect Control Hosp Epidemiol.* 28(6):684-8.

Infusion Therapy

Stopcocks for IV Therapy and Pressure Monitoring

BD Connecta™ - 3-way Stopcock

BD Connecta™ is the result of BD's continuous search for improvement and commitment to risk management in IV therapy.

BD Connecta™ offers the possibility to feel whether the stopcock is completely open or closed through the on/off snaps.

The new grade of polycarbonate and high standard moulding technique contribute to improved resistance to cracking.

- Flat clean surface
- Low priming volume (dead space)
- Straight two-way fluid path
- Material allows for use of lipids and chemotherapeutic agents
- Sterilization method is Radiation
- Sterile, single use
- CE marked



Description	Extension Tubing Length (cm)	Priming Volume (ml)	Shelf Box / Shipping Case	Reference Number
BD Connecta™ 3-way stopcock				
White (without color coding pegs)	–	0,22	100 / 500	394600
White (with color coding pegs red and blue)	–	0,22	100 / 500	394601
Blue	–	0,22	100 / 500	394602
Red	–	0,22	100 / 500	394605
BD Connecta™ 3-way stopcock with extension				
White	10	0,80	50 / 250	394995
Blue	10	0,80	50 / 250	394997
White	25	1,50	20 / 100	394926
White	50	2,80	20 / 100	394951
White	100	5,50	20 / 100	394961
BD Connecta™ 3-way stopcock with extension and valve				
White	7	0,80	50 / 250	394945
White	25	1,50	20 / 100	394936
White	100	5,50	20 / 100	394971
BD Connecta™ 3-way stopcock with low volume extension tube				
White	15	0,30	50 / 250	394982
White	30	0,45	50 / 250	394983

Infusion Therapy

Infusion Accessories

BD FloPro™ - Transfer Needle

The BD FloPro™ Transfer needle allows convenient, economic and safe mixing of pharmaceuticals prior to administration. The two spikes are siliconized for easier penetration

- Sterilization method is Ethylene Oxide
- Sterile, single use



Description	Shelf Box / Shipping Case	Reference Number
BD FloPro™ transfer needle	100 / 100	397015

Caps and Plugs

- Latex free caps for ease of application
- Designed to minimize touch-contamination risk
- Tray packaging reduces waste compared to single-sterile caps
- For Luer Lock and Luer Slip connections
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Connection	Color	Shelf Box / Shipping Case	Reference Number
Caps and Plugs				
BD Luer-Lok™ plug (tray of 4)	male	white	200 / 2000	394080
Combi plugs	male / female	orange	100 / 5000	394075
Rubber membrane Luer Lock plug	male	yellow	100 / 5000	394074
Protection cap for Luer Slip Syringes	male / female	yellow	100 / 5000	394090
BD Luer Tip Caps (tray of 10)	female	black	200 / 2000	302834

Infusion Therapy

One-way Valve

BD™ One-way Valve

When more than one IV line is connected to a single access point, back tracking of fluids can occur if one line has no flow or a slower flow of fluid running through it.

The fluid will back track as it will take the path of least resistance.

- This one-way Luer Lock valve is a mechanical device which allows fluid to flow through IV infusion lines in one way only

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Shelf Box / Shipping Case	Reference Number
BD™ One-way Valve	200 / 200	397019

Infusion Therapy

Skin Prep

BD Persist™ Plus

INS Standards (1998, S53) recommend antimicrobial solution containers be in a single-unit of use and that they should be discarded after individual use.

- Easy, one-step application
- Fast drying - 45 seconds vs. 3-5 minutes for povidone-iodine¹
- Unique green color
 - enhances visualisation of the zone of I.V. access
 - highlights the protective barrier
 - does not mask the vein
- Non-staining* and does not bond to uniforms, sheets or drapes

- Pleasantly scented, enhancing clinical and patient acceptance
- Available in aseptic, single use packing

¹ Maki DG, Knasinski V, Narans LL, Gordon BJ (2001) A randomized trial of a novel 1% chlorhexidine-75% alcohol tincture versus 10% povidone-iodine for cutaneous disinfection with vascular catheters. In: Proceedings and Abstracts 11th Annual Scientific Meeting, Society for Healthcare Epidemiology of America.

* when used correctly



Description	Shelf Box / Shipping Case	Reference Number
BD Persist™ Plus Swab Stick - Single	50 / 500	386411
BD Persist™ Plus Swab Stick - Triple Protective Barrier	25 / 250	386412

2011 CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections, Skin Preparation²

- Prepare clean skin with an antiseptic (70% alcohol, tincture of iodine, or alcoholic chlorhexidine gluconate solution) before peripheral venous catheter insertion category IB.
 - alcohol to prepare clean skin. Unresolved issue.
- Prepare skin with a > 0,5% chlorhexidine preparation with alcohol before central venous catheter and peripheral arterial catheter insertion and during dressing changes. If there is a contraindication to chlorhexidine, tincture of iodine, an iodophor, or 70 % alcohol can be used as alternatives. Category IA.
 - No recommendation can be made for the safety or efficacy of chlorhexidine in infants aged < 2 months. Unresolved issue.
 - Antiseptics should be allowed to dry according to the manufacturer's recommendation prior to placing the catheter. Category IB.
- No comparison has been made between using chlorhexidine preparation with alcohol and povidone-iodine in

² CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011

2012 JCI Guidelines - Preventing Central Line-Associated Bloodstream Infections

Recommendations for skin antiseptics prior to CVC insertion and during dressing changes³

- Apply antiseptics to clean skin
- Apply chlorhexidine/alcohol in a concentration greater than 0.5% in alcohol
 - If there is a contraindication to chlorhexidine, apply tincture of iodine, an iodophor, or alcohol as an alternative.
 - Allow the antiseptic solution to dry before placing the catheter.

³ JCI Guidelines 2012

Peripheral Vascular Access Systems

Closed IV Catheter Systems

BD Nexiva™ - Closed IV Catheter System

- BD Nexiva™ closed IV Catheter System is a pre-assembled configuration which combines separate IV components into a single device
- It promotes greater efficiency and convenience with less handling time, and reduces the potential for blood leakage, helping to minimize clean-up
- Pre-attached extension tubing and soft wings reduce manipulation at the insertion site to help decrease the risk of mechanical phlebitis and irritation to the vein

Improved first-stick success

- With its compact size and shape, the BD Nexiva™ Closed IV Catheter System accommodates your preferred insertion technique

BD Instaflash™ Needle Technology provides immediate visual confirmation of vessel entry along the cannula. These features increase the potential for a successful insertion

High Pressure Extension Set

- The pre-attached high pressure extension set is compatible with power injectors used with contrast-enhance CT scans*

- This preassembled system is designed to withstand pressures of up to a maximum of 300 PSI

BD Q-Syte™ Luer Access Split Septum

- A split septum needleless access system has 64%-70% lower CRBSI rates than mechanical valves^{1,2}
- Patients are three times more likely, on average, to develop a CRBSI with the use of a mechanical valve versus a split septum needleless access system^{1,2}

Built in Stabilisation Platform

- The built in stabilisation platform is soft and flexible, and is designed to help minimise cannula movement in the vessel which may increase cannula dwell time.

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- Manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



*24 gauge catheter system should not be used with power injectors

² Rupp ME, Sholtz LA, Jourdan DR et al. (2007): Outbreak of bloodstream infection temporally associated with the use of an intravascular needleless valve. *Clin Infect Dis.* 44(11):1408-14.

³ Salgado CD, Chinnes L, Paczesny TH, Cantey JR (2007): Increased rate of catheter-related bloodstream infection associated with use of a needleless mechanical valve device at a long-term acute care hospital. *Infect Control Hosp Epidemiol.* 28(6):684-8.


BD Instaflash™ Technology - Designed to increase first-stick success

A successful first-stick is important for both the user and the patient, especially when accessing small and compromised veins. BD Instaflash™ needle technology is designed to reduce hit-and-miss insertion by immediately confirming vessel entry through a unique notch in the needle shaft.

Additional assurance of vessel entry is provided through continuous flashback and clear visibility of blood flow along the extension tubing.

Passive needle protection device reduces the potential for needlestick injuries. The needle tip is automatically shielded when withdrawn from the cannula to ensure compliance. The shield remains in place and cannot be bypassed.

Closed fluid flow path reduces the potential for blood leakage, thereby reducing clinicians' exposure to bloodborne pathogens.



BD Instaflash™ Technology. Designed to reduce hit-and-miss insertion.

① BD Instaflash™ Technology



Peripheral Vascular Access Systems

Closed IV Catheter Systems

BD Nexiva™ - Closed IV Catheter System

Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate H ₂ O (ml/min)	Priming Volume (ml)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Nexiva™ Single Port							
24 G	0,7	14	16		yellow	20 / 80	383510
24 G	0,7	19	15		yellow	20 / 80	383511
22 G	0,9	25	30		blue	20 / 80	383512
20 G	1,1	25	54		pink	20 / 80	383516
20 G	1,1	32	52		pink	20 / 80	383517
20 G	1,1	45	48		pink	20 / 80	383518
18 G	1,3	32	81		green	20 / 80	383519
18 G	1,3	45	76		green	20 / 80	383520
BD Nexiva™ Dual Port "Y" with one BD Q-Syte™							
24 G	0,7	14	19	0,4	yellow	20 / 80	383640
24 G	0,7	19	18	0,4	yellow	20 / 80	383641
22 G	0,9	25	33	0,4	blue	20 / 80	383642
20 G	1,1	25	61	0,5	pink	20 / 80	383646
20 G	1,1	32	58	0,5	pink	20 / 80	383647
20 G	1,1	45	51	0,5	pink	20 / 80	383648
18 G	1,3	32	84	0,5	green	20 / 80	383649
18 G	1,3	45	79	0,6	green	20 / 80	383650
BD Nexiva™ Dual Port "Y" with two BD Q-Syte™							
24 G	0,7	14	19	0,4	yellow	20 / 80	383530
24 G	0,7	19	18	0,4	yellow	20 / 80	383531
22 G	0,9	25	33	0,4	blue	20 / 80	383532
20 G	1,1	25	61	0,5	pink	20 / 80	383536
20 G	1,1	32	58	0,5	pink	20 / 80	383537
18 G	1,3	32	84	0,5	green	20 / 80	383539
18 G	1,3	45	79	0,6	green	20 / 80	383540
BD Nexiva™ with BD Connecta™ Stopcock							
24 G	0,7	19	18	0,6	yellow	20 / 80	383681
22 G	0,9	25	33	0,6	blue	20 / 80	383682
20 G	1,1	25	61	0,8	pink	20 / 80	383686
20 G	1,1	32	58	0,8	pink	20 / 80	383687
18 G	1,3	32	84	0,8	green	20 / 80	383689
BD Nexiva™ with BD Connecta™ Stopcock and BD Q-Syte™							
24 G	0,7	14	19	0,6	yellow	20 / 80	383660
24 G	0,7	19	18	0,6	yellow	20 / 80	383661
22 G	0,9	25	33	0,6	blue	20 / 80	383662
20 G	1,1	25	61	0,8	pink	20 / 80	383666
20 G	1,1	32	58	0,8	pink	20 / 80	383667
20 G	1,1	45	51	0,8	pink	20 / 80	383668
18 G	1,3	32	84	0,8	green	20 / 80	383669
18 G	1,3	45	79	0,8	green	20 / 80	383670

Peripheral Vascular Access Systems


Closed IV Catheter Systems

BD Nexiva™ Diffusics™ - Closed IV Catheter System Radiographic Power Injection

BD Nexiva™ Diffusics™ is an innovative addition to the already proven and effective IV catheter system, BD Nexiva™. It is a pre-assembled closed IV catheter system which minimizes blood exposure, improves catheter stabilization, enhances insertion success with BD Instaflash™ needle technology whilst being specifically designed for radiology needs.

BD Nexiva™ Diffusics™ was specifically engineered to deliver on today's power injection and infusion therapy needs. It is built to allow clinicians easy IV access and to deliver the desired injection protocol flow rates in diagnostic imaging.

BD Nexiva™ Diffusics™ has three laser-cut teardrop holes optimally shaped and located in the catheter tip to maximize the diffusion of contrast into the vein. The BD Diffuser tip acts to reduce the intensity




Increased diffusion of contrast due to teardrop shaped holes in the catheter tip

of flow exiting the catheter tip, which may decrease stress on the vein wall.

This all-in-one BD IV access system includes the IV catheter, extension set and stabilization platform. It is compatible with power injector's 325 psi setting **AND** can meet all high flow rate injection protocols.

- Now even smaller IV catheters gauge sizes (22 and 24G) can meet the needs of faster protocols
- Reduces the forces that can cause catheter motion in the vein¹
- Significantly less impact to the vein during power injection²
- Helps protect clinicians from blood exposure during IV cannulation
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



BD Instaflash™ Technology

Gauge Size	Extension Tube I.D. (mm)	Catheter Length (mm)	max. CT Flow Rate (ml/sec)*	Gravity Flow Rate (ml/h)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Nexiva™ Diffusics™ Closed IV Catheter System for Radiographic Power Injection							
24 G	1,22	19	3.0	1235	yellow	20 / 80	383691
22 G	1,65	25	6.5	2690	blue	20 / 80	383692
20 G	1,65	25	10.0	4050	pink	20 / 80	383693
20 G	1,65	32	10.0	3850	pink	20 / 80	383694
18 G	1,65	32	15.0	5400	green	20 / 80	383695

* see instructions for use for all power injection guidelines, cautions and warnings



Scan this code with your phone or tablet
www.bd.com/diffusics



¹ An in-vitro study measuring IV catheter recoil forces during power injection of iodinated contrast media. 2011 BD ITS0164-1B (11/11).

² An in-vitro assessment of diffuser efficiency during power injection of iodinated contrast media using the BD Nexiva™ Diffusics™ Closed IV Catheter System. 2011 BD ITS0165-1B (11/11).

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Venflon™ Pro Safety - Peripheral Safety IV Catheter with Injection Valve

BD Venflon™ Pro Safety needle protected IV Cannula combines all the advantages of the original BD Venflon™ Pro, with second generation safety protection features against needlestick injuries.

- Familiar ported design with 'snap-cap' aids non-touch technique which helps prevent port being inadvertently left open
- Sharp needle technology provides both ease of penetration and improved patient comfort
- Flexible wings stabilize the cannula in the vein and provide comfortable securement for the patient
- BD Venflon™ Pro Safety is designed for use in all departments and provides a full range of features to suit most clinical indications
- Grip is ergonomically designed to allow for a variety of user techniques and to aid insertion

- Needle shield
 - once activated, the needle tip is fully encapsulated inside the protection mechanism, which is designed to minimise the risk of scratching yourself
 - the unique design captures residual blood safely within the device, thereby minimizing the risk of unanticipated blood splatter during and after activation
 - with needle safely encased in smooth edged plastic shield
- Sterilization method is Radiation
- Sterile, single use
- CE marked
- Manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Venflon™ Pro Safety						
22 G	0,9	25	42	blue	50 / 500	393222
20 G	1,1	32	67	pink	50 / 500	393224
18 G	1,3	32	103	green	50 / 500	393226
18 G	1,3	45	103	green	50 / 500	393227
17 G	1,5	45	133	white	50 / 500	393228
16 G	1,8	45	236	grey	50 / 500	393229
14 G	2,0	45	270	orange	50 / 500	393230



reddot design award
best of the best 2008

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Insyte™ Autoguard™ / BD Insyte-N™ Autoguard™ - Peripheral Safety IV Catheter


BD knows the importance to both clinician and patient of achieving first-stick success which is why we developed the BD Insyte™ Autoguard™ shielded IV catheter.

BD Insyte™ Autoguard™ shielded IV catheters are designed to give you the satisfaction that comes from getting a successful first stick.

- Patented safety mechanism activation
- The push-button holds the needle unit in the extended position until the user activates the mechanism
- Pressing the push-button releases the spring so that the needle/flashback chamber retracts into the safety barrel.
- Luer Lock compatible

BD Insyte-N™ Autoguard™

- 14 mm long, the shortest on the market to accommodate small vein segments and provides a greater choice of insertion sites
- 24G for catheterisation in every situation encountered in neonatology, paediatrics, oncology and elderly patient
- Shorter length provides increased flow rates
- Available with or without wings
- Recommended uses:
 - fragile veins, e.g. neonatology
 - low venous pressure
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



BD Instafash™ Technology

Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Insyte-N™ Autoguard™ winged						
24 G*	0,7	14	20	yellow	50 / 200	381911
BD Insyte™ Autoguard™ winged						
24 G*	0,7	19	20	yellow	50 / 200	381912
22 G*	0,9	25	35	blue	50 / 200	381923
20 G*	1,1	25	65	pink	50 / 200	381933
20 G*	1,1	30	60	pink	50 / 200	381934
20 G*	1,1	48	55	pink	50 / 200	381937
18 G	1,3	30	105	green	50 / 200	381944
18 G	1,3	48	95	green	50 / 200	381947
16 G	1,7	30	220	grey	50 / 200	381954
16 G	1,7	45	205	grey	50 / 200	381957
BD Insyte-N™ Autoguard™ without wings						
24 G*	0,7	14	20	yellow	50 / 200	381811
BD Insyte™ Autoguard™ without wings						
24 G*	0,7	19	20	yellow	50 / 200	381812
22 G*	0,9	25	35	blue	50 / 200	381823
20 G*	1,1	25	65	pink	50 / 200	381833
20 G*	1,1	30	60	pink	50 / 200	381834
20 G*	1,1	48	55	pink	50 / 200	381837
18 G	1,3	30	105	green	50 / 200	381844
18 G	1,3	48	95	green	50 / 200	381847
16 G	1,7	30	220	grey	50 / 200	381854
16 G	1,7	45	205	grey	50 / 200	381857
14 G	2,1	45	330	orange	50 / 200	381867

* with BD Instafash™ Technology

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Insyte™ Autoguard™ BC - with Blood Control Technology

BD Insyte™ Autoguard™ BC is proven to reduce the risk of blood exposure by 95%.¹

Reduce cross contamination risks

- Because there is dramatically less blood leakage, BD Insyte™ Autoguard™ BC can mean less potential cross contamination
- No spills means there's no blood to get on gloves, clothing, shoes, bedding, floors, or equipment

Safe and clinically proven

- New blood-control technology proven to reduce the risk of blood exposure by 95%¹
- Push-button shielding technology proven to reduce needlestick injuries by 95%²
- Greater first-stick success for difficult or compromised veins³

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- BD Vialon™ Biomaterial—softens up to 70% in the vessel, resulting in longer dwell times. Reduces the chance of mechanical phlebitis by up to 50%^{4, 5, 6, 7}

¹ Data on file

² Mendelson MH, Lin-Chen BY, Finkelstein-Blond L, Bailey E, Kogan G. Evaluation of a safety IV catheter (IVC) (Becton Dickinson, Insyte Autoguard): final report [abstract]. Eleventh Annual Scientific Meeting Society for Healthcare Epidemiology of America. 2001.

³ Bausone-Gazda D, Lefaiver C, Walters S. A randomized controlled trial to compare the complications of 2 peripheral intravenous catheter-stabilization systems. *J Infus Nurs.* 2010;33(6):371-384.

⁴ Maki DG, Ringer M. Risk factors for infusion-related phlebitis with small peripheral venous catheters. *Annals of Internal Medicine.* 1991;114:845-854.


⁵ Gaukroger PB, Roberts JG, Manners TA. Infusion thrombophlebitis: a prospective comparison of 645 Vialon and Teflon cannulae in anaesthetic and postoperative use. *Anaesthesia and Intensive Care.* 1998;16:265-271.

⁶ Stanley MD, Meister E, Fuschuber K. Infiltration during intravenous therapy in neonates: comparison of Teflon® and Vialon® catheters. *Southern Medical Journal.* 1992;85:883-886.

⁷ McKee JM, Shell JA, Warren TA, Campbell VP. Complications of intravenous therapy: a randomized prospective study — Vialon vs. Teflon. *J Infus Nurs.* 1989;12:288-295.



BD Blood Control Technology – stops blood so it doesn't escape the hub after insertion



BD Instaflash™ Technology

Gauge Size	Outer diameter (mm)	Inner Diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Insyte™ Autoguard™ winged BC							
24 G*	0,66 - 0,74	0,48 - 0,56	19	20	yellow	50 / 200	382912
22 G*	0,84 - 0,91	0,61 - 0,69	25	37	blue	50 / 200	382923
20 G*	1,04 - 1,12	0,76 - 0,84	25	63	pink	50 / 200	382933
20 G*	1,04 - 1,12	0,76 - 0,84	30	61	pink	50 / 200	382934
20 G*	1,04 - 1,12	0,76 - 0,84	48	54	pink	50 / 200	382937
18 G	1,27 - 1,35	0,94 - 1,02	30	95	green	50 / 200	382944
18 G	1,27 - 1,35	0,94 - 1,02	48	87	green	50 / 200	382947
16 G	1,70 - 1,78	1,32 - 1,40	30	193	grey	50 / 200	382954
16 G	1,70 - 1,78	1,32 - 1,40	45	185	grey	50 / 200	382957
BD Insyte™ Autoguard™ BC without wings							
24 G*	0,66 - 0,74	0,48 - 0,56	19	20	yellow	50 / 200	381012
22 G*	0,84 - 0,91	0,61 - 0,69	25	37	blue	50 / 200	381023
20 G*	1,04 - 1,12	0,76 - 0,84	25	63	pink	50 / 200	381033
20 G*	1,04 - 1,12	0,76 - 0,84	30	61	pink	50 / 200	381034
20 G*	1,04 - 1,12	0,76 - 0,84	48	54	pink	50 / 200	381037
18 G	1,27 - 1,35	0,94 - 1,02	30	95	green	50 / 200	381044
18 G	1,27 - 1,35	0,94 - 1,02	48	87	green	50 / 200	381047
16 G	1,70 - 1,78	1,32 - 1,40	30	193	grey	50 / 200	381054
16 G	1,70 - 1,78	1,32 - 1,40	45	185	grey	50 / 200	381057

* with BD Instaflash™ Technology

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Saf-T-Intima™ – Integrated Safety IV Catheter System

- BD Saf-T-Intima™ is equipped with a needle protection device that withdraws the needle into a safety barrel as soon as the needle is removed. This minimises the risk of accidental needle stick injuries
- BD Saf-T-Intima™ is fitted with a BD Luer-Lok™ connector for connection to any piece of equipment fitted with a Luer Adapter
- The BD Saf-T-Intima™ “Y” type adapter allows that the injection site can be used for bolus administration of medication
- Removable injection port
- Blood contact is avoided during insertion
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- Manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



Gauge Size	Outer Diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Subcutaneous Application	Color Code	Shelf Box / Shipping Case	Reference Number
BD Saf-T-Intima™ with Y-adapter							
24 G	0,7	19	22	✓	yellow	25 / 200	383319
22 G	0,9	19	36	✓	blue	25 / 200	383329
20 G	1,1	25	55	–	pink	25 / 200	383339
18 G	1,3	25	82	–	green	25 / 200	383348
BD Saf-T-Intima™ without Y-adapter							
24 G	0,7	19	22	✓	yellow	25 / 200	383318
22 G	0,9	19	36	✓	blue	25 / 200	383328
20 G	1,1	25	55	–	pink	25 / 200	383338

BD Saf-T-Intima™ is now approved by the British Standard Institution (BSI) for subcutaneous use, it provides patient satisfaction, HCW safety and cost efficiency. The use of the BD Saf-T-Intima™ Integrated Safety IV Catheter System for subcutaneous use is an extension to this product's indications for use.

Subcutaneous use was not within the scope of the original BD Saf-T-Intima™ CE Certificate. This extension required full assessment by BD's notified body before BD Saf-T-Intima™ could be used subcutaneously within the European Union. BD has sought to extend the intended use of BD Saf-T-Intima™ because there is a

growing practice of subcutaneous infusion, as it can potentially offer a better, more comfortable solution for treating rehydration, pain relief and end of life care, than the traditional IV route.



BD Saf-T-Intima™ Integrated Safety IV Catheter System is approved for subcutaneous infusion therapies by BSI.

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Venflon™ Pro - Peripheral IV Catheter with Injection Valve

- The improved design of the injection port cap ensures easy opening and secure closing, whilst avoiding touch contamination of the injection port and valve
- Improved insertion characteristics as BD Vialon™ make it possible to manufacture a catheter with a better tapered, stronger and more reliable tip geometry which improves the chances of a first cannulation attempt success and also reduces the risk of thrombophlebitis
- Although the catheter material is rigid before insertion, the qualities of BD Vialon™ allow for softening of the catheter material once it is inside the body, again reducing the risk of thrombophlebitis
- There is an improved flow rate when compared with FEP polymer catheters as BD Vialon™ allows for a larger internal diameter
- BD Luer-Lok™ compatible
- Sterilization method is Radiation
- Sterile, single use
- CE marked
- Manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Venflon™ Pro						
22 G	0,9	25	42	blue	50 / 500	393202
20 G	1,1	32	67	pink	50 / 500	393204
18 G	1,3	32	103	green	50 / 500	393206
18 G	1,3	45	103	green	50 / 500	393207
17 G	1,5	45	133	white	50 / 500	393208
16 G	1,8	45	236	grey	50 / 500	393209
14 G	2,0	45	270	orange	50 / 500	393210

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Venflon™ - Peripheral IV Catheter with Injection Valve

- BD Venflon™ intravenous catheter with integrated injection port and valve for peripheral catheterisation
- Permits needle-free administration of bolus medications during infusion
- Back-cut needle bevel and tapered end of the PTFE catheter allows easier insertion and minimizes the risks of vein damage
- Wing position reduces the risk of bacterial contamination by preventing any direct contact between the rear of the catheter hub and the skin. Utilizing the wings, immobility of the catheter can be further ensured by the application of the BD Veca-C™ and BD Vecafix™ IV dressing
- BD Luer-Lok™ compatible
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Venflon™						
22 G	0,8	25	31	blue	50 / 500	391451
20 G	1,0	32	54	pink	50 / 500	391452
18 G	1,2	32	80	green	50 / 500	391457
18 G	1,2	45	80	green	50 / 500	391453
17 G	1,4	45	125	white	50 / 500	391454
16 G	1,7	45	180	grey	50 / 500	391455
14 G	2,0	45	270	orange	50 / 500	391456

Peripheral Vascular Access Systems

Indwelling Cannulae


BD Insyte™ and BD Insyte-N™ - Peripheral IV Catheter with and without Wings

BD Insyte™

- Easy insertion
 - the strength of BD Vialon™ makes it possible to manufacture a catheter with a finer, stronger, more reliable tip, giving it excellent penetration characteristics
- The catheter is rigid before insertion but becomes soft inside the body. This unique feature adds to the biocompatibility of BD Vialon™ and results in a lower risk of complications
- Flow rate: BD Insyte™ is a thin-walled catheter with a higher flow rate than other catheters of the same gauge. BD Vialon™ has good kink resistance and kink recovery to ensure an optimum flow rate
- Recommended uses:
 - Oncology, low venous pressure, fragile veins

BD Insyte-N™

- 14 mm long, the shortest on the market to accommodate small vein segments and provides a greater choice of insertion sites
- 24G for catheterization in every situation encountered in neonatology, paediatrics, oncology and geriatrics.
- BD Luer-Lok™ compatible
- Recommended uses:
 - fragile veins, e.g. neonatology
 - low venous pressure
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- Manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



BD Instaflash™ Technology

Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Insyte-N™ with wings						
24 G*	0,7	14	25	yellow	50 / 200	381311
BD Insyte-W™ with wings						
24 G*	0,7	19	25	yellow	50 / 200	381312
22 G*	0,9	25	35	blue	50 / 200	381323
20 G*	1,1	30	60	pink	50 / 200	381334
20 G*	1,1	48	55	pink	50 / 200	381337
18 G	1,3	30	105	green	50 / 200	381344
18 G	1,3	48	95	green	50 / 200	381347
16 G	1,7	45	205	grey	50 / 200	381357
BD Insyte-N™ without wings						
24 G*	0,7	14	25	yellow	50 / 200	381211
BD Insyte™ without wings						
24 G*	0,7	19	25	yellow	50 / 200	381212
22 G*	0,9	25	35	blue	50 / 200	381223
20 G*	1,1	30	60	pink	50 / 200	381234
20 G*	1,1	48	55	pink	50 / 200	381237
18 G	1,3	30	105	green	50 / 200	381244
18 G	1,3	48	95	green	50 / 200	381247
16 G	1,7	45	205	grey	50 / 200	381257
14 G	2,1	45	330	orange	50 / 200	381267

* with BD Instaflash™ Technology

Peripheral Vascular Access Systems

Indwelling Cannulae

BD Angiocath™ - Pheripheral IV Catheter

- Catheter made of radiopaque FEP Polymer.
- Catheter with thermo-moulded bevel ensures excellent penetration properties
- Triple-bevel needle
- Transparent clear flashback chamber
- Complete gauge range from 10G x 76 mm to 24G x 19 mm
- Available in longer lengths (14G and 16G available with 83 mm and 133 mm long catheters)
- BD Luer-Lok™ compatible
- Recommended uses
 - Infusion
 - Transfusion
 - Angiography
- Sterilization method is Ethylene Oxide
- CE marked
- Sterile, single use
- NB: Diameters 12G and 10G permit very large volume emergency infusions



Gauge Size	Outer Diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min.)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Angiocath™						
24 G	0,7	19	17	yellow	50 / 200	381112
22 G	0,9	25	28	blue	50 / 200	381123
20 G	1,1	30	49	pink	50 / 200	381134
20 G	1,1	48	42	pink	50 / 200	381137
18 G	1,3	30	86	green	50 / 200	381144
18 G	1,3	48	79	green	50 / 200	381147
16 G	1,7	30	156	grey	50 / 200	381154
16 G	1,7	48	147	grey	50 / 200	381157
16 G	1,7	83	127	grey	10 / 50	382258
16 G	1,7	133	108	grey	10 / 50	382259
14 G	2,1	30	295	orange	50 / 200	381164
14 G	2,1	48	277	orange	50 / 200	381167
14 G	2,1	83	249	orange	10 / 50	382268
14 G	2,1	133	219	orange	10 / 50	382269
12 G	2,8	76	449	light blue	10 / 50	382277
10 G	3,4	76	609	tan	10 / 50	382287

Peripheral Vascular Access Systems

Indwelling Cannulae for Neonatal and Pediatric Use / IV Dressing

BD Neoflon™

- 24G and 26G PTFE catheter
- BD Neoflon™ is specially designed for the puncture of small fragile vessels, frequently encountered in paediatrics or elderly patients
- The unique design of the detachable wing holder makes the catheter easier to hold, using either a ported or non-ported grip, despite its small size
- The flexible wings facilitate secure fixation, even on the hand or arm of the neonate without the need to immobilize the limb whilst at the same time reducing the risk of mechanical phlebitis
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Outer Diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Neoflon™ with wings						
24 G	0,7	19	13	yellow	50 / 500	391350
26 G	0,6	19	13	violet	50 / 500	391349

BD Veca-C™ and BD Vecafix™ - IV Dressing

BD Veca-C™ and BD Vecafix™ are sterile IV dressings for optimum fixation of ported and non-ported IV catheters.

- BD Veca-C™ and BD Vecafix™ protects the insertion site and allows the site to breathe
- The wound pad absorbs exudate at the insertion site, keeping the site dry
- A comfort pad, to be placed under the wings of the IV catheter, is included in the pack
- In addition, BD Veca-C™ has a transparent window of polyurethane that allows for visualization of phlebitis along the catheter track
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- The packing of the BD Veca-C™ dressing contains small amounts of latex. The product itself however, is latex free



Description	Dimension (cm)	Shelf Box / Shipping Case	Reference Number
BD Veca-C™ IV dressing with transparent window			
	7,5 x 6,0	50 / 600	392020
BD Vecafix™ IV dressing without transparent window			
	7,5 x 6,0	50 / 600	392015

Peripheral Vascular Access Systems

Obturator

Obturator

- The Obturator protects the inner lumen of the catheter from coagulation and thrombus formation after discontinuation of infusion
- The Obturator is color coded to match the corresponding catheter gauge
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked

Alternatively, we recommend the flushing of the catheter with prefilled flushing devices in order to keep the catheter open.

Please refer to BD PosiFlush™.

Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Indwelling Cannula		Obturator	
					Shelf Box / Shipping Case	Reference Number	Shelf Box / Shipping Case	Reference Number
BD Venflon™								
22 G	0,8	25	31	blue	50 / 500	391451	50 / 500	394251
20 G	1,0	32	54	pink	50 / 500	391452	50 / 500	394252
18 G	1,2	32	80	green	50 / 500	391457	50 / 500	394257
18 G	1,2	45	80	green	50 / 500	391453	50 / 500	394253
17 G	1,4	45	125	white	50 / 500	391454	50 / 500	394254
16 G	1,7	45	180	grey	50 / 500	391455	50 / 500	394255
14 G	2,0	45	270	orange	50 / 500	391456	50 / 500	394256
BD Venflon™ Pro								
22 G	0,9	25	42	blue	50 / 500	393202	50 / 500	394251
20 G	1,1	32	67	pink	50 / 500	393204	50 / 500	394252
18 G	1,3	32	103	green	50 / 500	393206	50 / 500	394257
18 G	1,3	45	103	green	50 / 500	393207	50 / 500	394253
17 G	1,5	45	133	white	50 / 500	393208	50 / 500	394254
16 G	1,8	45	236	grey	50 / 500	393209	50 / 500	394255
14 G	2,0	45	270	orange	50 / 500	393210	50 / 500	394256
BD Venflon™ Pro Safety								
22 G	0,9	25	42	blue	50 / 500	393222	50 / 500	394251
20 G	1,1	32	67	pink	50 / 500	393224	50 / 500	394252
18 G	1,3	32	103	green	50 / 500	393226	50 / 500	394257
18 G	1,3	45	103	green	50 / 500	393227	50 / 500	394253
17 G	1,5	45	133	white	50 / 500	393228	50 / 500	394254
16 G	1,8	45	236	grey	50 / 500	393229	50 / 500	394255
14 G	2,0	45	270	orange	50 / 500	393230	50 / 500	394256

Peripheral Vascular Access Systems

Obturator

Gauge Size	Outer diameter (mm)	Catheter Length (mm)	Flow Rate (ml/min)	Color Code	Indwelling Cannula		Obturator	
					Shelf Box / Shipping Case	Reference Number	Shelf Box / Shipping Case	Reference Number
BD Insyte-W™ with wings								
24 G*	0,7	19	25	yellow	50 / 200	381312	50 / 500	394287
22 G*	0,9	25	35	blue	50 / 200	381323	50 / 500	394280
20 G*	1,1	30	60	pink	50 / 200	381334	50 / 500	394281
20 G*	1,1	48	55	pink	50 / 200	381337	50 / 500	394282
18 G	1,3	30	105	green	50 / 200	381344	50 / 500	394283
18 G	1,3	48	95	green	50 / 200	381347	50 / 500	394284
16 G	1,7	45	205	grey	50 / 200	381357	50 / 500	394285
BD Insyte™ without Wings								
24 G*	0,7	19	25	yellow	50 / 200	381212	50 / 500	394287
22 G*	0,9	25	35	blue	50 / 200	381223	50 / 500	394280
20 G*	1,1	30	60	pink	50 / 200	381234	50 / 500	394281
20 G*	1,1	48	55	pink	50 / 200	381237	50 / 500	394282
18 G	1,3	30	105	green	50 / 200	381244	50 / 500	394283
18 G	1,3	48	95	green	50 / 200	381247	50 / 500	394284
16 G	1,7	45	205	grey	50 / 200	381257	50 / 500	394285
14 G	2,1	45	330	orange	50 / 200	381267	50 / 500	394286
BD Insyte™ Autoguard™ with wings								
24 G*	0,7	19	20	yellow	50 / 200	381912	50 / 500	394287
22 G*	0,9	25	35	blue	50 / 200	381923	50 / 500	394280
20 G*	1,1	30	60	pink	50 / 200	381934	50 / 500	394281
20 G*	1,1	48	55	pink	50 / 200	381937	50 / 500	394282
18 G	1,3	30	105	green	50 / 200	381944	50 / 500	394283
18 G	1,3	48	95	green	50 / 200	381947	50 / 500	394284
16 G	1,7	45	205	grey	50 / 200	381957	50 / 500	394285
BD Insyte™ Autoguard™ without wings								
24 G*	0,7	19	20	yellow	50 / 200	381812	50 / 500	394287
22 G*	0,9	25	35	blue	50 / 200	381823	50 / 500	394280
20 G*	1,1	30	60	pink	50 / 200	381834	50 / 500	394281
20 G*	1,1	48	55	pink	50 / 200	381837	50 / 500	394282
18 G	1,3	30	105	green	50 / 200	381844	50 / 500	394283
18 G	1,3	48	95	green	50 / 200	381847	50 / 500	394284
16 G	1,7	45	205	grey	50 / 200	381857	50 / 500	394285
14 G	2,1	45	330	orange	50 / 200	381867	50 / 500	394286

* with BD Instaflash™ Technology

Peripheral Vascular Access Systems

BD Vialon™ Biomaterial

BD Vialon

BD Vialon™ Biomaterial is clinically proven to optimize patient outcomes

BD Vialon™ biomaterial is a unique, proprietary biomaterial, developed especially for vascular access.

Multiple clinical investigations conducted over ten years have proven that IV catheters made of BD Vialon™ biomaterial are easier to insert, while providing longer dwell times with improved patient outcomes.^{1,2,3,4}

The insertion characteristics of BD Vialon™ biomaterial catheters are as unique as the material itself. Superior tip integrity, an ultra-smooth micro-surface and an exclusive lubrication process enable the easiest insertion possible, making BD Vialon™ biomaterial catheters easier to place and advance on a wider cross-section of patients.²

Softens inside the vessel

Another distinctive characteristic of BD Vialon™ biomaterial is its ability to soften up to 70% while in the vessel.⁵

By conforming to the shape of the vein, catheters made of BD Vialon™ biomaterial reduce the chance of irritation to the vessel wall, which in turn reduces incidents of mechanical phlebitis, a major complication of IV therapy.

Optimum reduced risk of mechanical phlebitis and IV infiltration

Vialon catheters changed every 3 or 4 days produce a lower incidence of phlebitis than would other FEP catheters changed every 2 days.¹

The reduction in risk of phlebitis was nearly 30% overall; in severe phlebitis, the reduction was nearly 50%.¹ Vialon cannulae are associated with a 46% lower incidence and severity of infusion thrombophlebitis.²

¹ Maki D, Ringer M. Risk Factors for Infusion-related Phlebitis with Small Peripheral Venous Catheters. *Annals of Internal Medicine*. 1991;114:845-854.

² Mendelson MH, Lin-Chen BY, Finkelstein-Blond L, Bailey E, Kogan G (2001): Evaluation of a safety IV catheter (IVC) (Becton Dickinson, Insyte Autoguard): In: Final report [abstract] 11th Annual Scientific Meeting, Society for Healthcare Epidemiology of America.


³ Bausone-Gazda D, Lefaiver CA, Walters SA (2010): A randomized controlled trial to compare the complications of 2 peripheral intravenous catheter-stabilization systems. *J Infus Nurs*. 33(6):371-84.

⁴ Maki DG, Ringer M (1991): Risk factors for infusion-related phlebitis with small peripheral venous catheters. A randomized controlled trial. *Ann Intern Med*. 114(10):845-54.


⁵ Gaukroger PB, Roberts JG, Manners TA (1988): Infusion thrombophlebitis: a prospective comparison of 645 Vialon and Teflon cannulae in anaesthetic and postoperative use. *Anaesth Intensive Care*. 16(3):265-71.

⁶ Stanley MD, Meister E, Fuschuber K (1992): Infiltration during intravenous therapy in neonates: comparison of Teflon and Vialon catheters. *South Med J*. 85(9):883-6

⁷ McKee JM, Shell JA, Warren TA, Campbell VP (1989): Complications of intravenous therapy: a randomized prospective study - Vialon vs. Teflon. *J Intraven Nurs*. 12(5):288-95.



The BD Vialon™ Optimum Kink Resistance BD Vialon biomaterial is markedly more flexible than other IV catheter materials.⁵ According to one study, "Kinking was not found in any of the Vialon cannulae after removal."



BD Vialon™

Products manufactured with BD Vialon™

BD Insyte™
BD Insyte-N™ Autoguard™
BD Insyte™ Autoguard™
BD Insyte™ Autoguard™ BC
BD Insyte-W™
BD Insyte-N™
BD Nexiva™
BD Nexiva™ Diffusics™
BD Saf-T-Intima™
BD Venflon™ Pro
BD Venflon™ Pro Safety

Peripheral Vascular Access Systems

BD Cannulae / Catheter Conversions

Color coding and corresponding sizes of Cannulae / Catheter


Outside diameter of catheter in mm	Range of actual outside diameter in mm	Color ^{1) 2)}	Gauge ³⁾
0,6	0,550 to 0,649	purple	26
0,7	0,650 to 0,749	yellow	24
0,8; 0,9	0,750 to 0,949	blue	22
1,0; 1,1	0,950 to 1,149	pink	20
1,2; 1,3	1,150 to 1,349	green	18
1,4; 1,5	1,350 to 1,549	white	17
1,6; 1,7; 1,8	1,550 to 1,849	grey	16
1,9; 2,0; 2,1; 2,2	1,850 to 2,249	orange	14
2,3; 2,4; 2,5	2,250 to 2,549	red	13
2,6; 2,7; 2,8	2,550 to 2,849	light blue	12
3,3; 3,4	3,250 to 3,549	tan	10

¹⁾ The color may be opaque or translucent.

²⁾ The color coding is usually applied to the catheter hub or to an integral fitting.


³⁾ The use of gauge numbers is optional.

Cannulae / Catheter Conversions Outer Diameters / Gauge-French Sizes



The definition of gauge measurement comes from the AAMI (Association for the Advancement of Medical Instrumentation - USA). This states how many cannulae fit into a tube with an inner diameter of 1" (1 inch = 25.4 mm).

The AAMI uses only even numbers (18, 20, 22, etc.) as measuring units. A similar measurement definition exists in Great Britain and this is known as the SWG (Standard Wire Gauge). The SWG uses consecutive numbers, i.e. 13 to 24 and is the more common method of measurement for catheters in Europe. The Charriere, commonly known as French, relates directly to the catheter size, i.e. 1 Ch = 0.33 mm. In the following table the measurements of the AAMI - SWG - French are shown in relation to the outer \varnothing in mm.



Flushing of Vascular Access Devices

BD PosiFlush™ - Prefilled Flushing Syringes [NaCl 0,9 %]

- Ready-to-use, terminally sterile prefilled flushing devices containing Sodium Chloride (isotonic NaCl 0.9%)
- Specially designed for flushing in-situ vascular access devices
- Special design of syringe barrel. The 3 and 5 ml syringes have the same barrel diameter as the 10 ml syringe in order to assure a standardized pressure during the flushing process
- BD Luer-Lok™ tip for secure connection to the catheter
- BD PosiFlush™ syringes are designed to eliminate syringe-induced blood reflux enhancing catheter maintenance protocols
- Size and shape of the tip cap allow easy handling, and help avoiding touch contamination of the syringe Luer connection
- Eliminates the use of needles, thus reducing the potential for accidental needle stick injuries
- Does not contain natural rubber latex
- Terminally steam sterilized - sterility assurance level of 10⁻⁶
- CE marked



*Average reflux as measured in 4Fr PICC; data on file at BD

Description	Shelf Box / Shipping Case	Reference Number
BD PosiFlush™ XS - externally sterile, for use in critical care areas for sterile procedures and within a sterile field		
Prefilled with 3 ml Saline Solution (NaCl 0,9 %)	30 / 240	306580
Prefilled with 5 ml Saline solution (NaCl 0,9 %)	30 / 240	306581
Prefilled with 10 ml Saline solution (NaCl 0,9 %)	30 / 240	306582
BD PosiFlush™ SP - Sterile Fluid Path, for general IV flush procedures		
Prefilled with 3 ml Saline Solution (NaCl 0,9 %)	30 / 480	306583
Prefilled with 5 ml Saline Solution (NaCl 0,9 %)	30 / 480	306584
Prefilled with 10 ml Saline solution (NaCl 0,9 %)	30 / 480	306585


Flushing of Vascular Access Devices

Consistent 10 ml diameter designed to lower the risk of catheter damage

Syringe size has an impact on the risk of catheter damage. Smaller diameter syringes generate greater amounts of pressure than larger diameter syringes.¹


Choose any size BD PosiFlush™ Syringe for lower pressure due to 10 ml syringe diameter (as compared to standard syringes).

All BD PosiFlush™ Syringes assure compliance with PICC manufacturer recommendations for flushing with 10 ml diameter syringes.



¹ Hadaway L (2006): Technology of flushing vascular access devices. J Infus Nurs. 29(3):137-45


BD PosiFlush™ Syringes are available in 3 sizes, all with standard 10 ml syringe diameter ...



Arterial Catheters

BD Insyte-A™

- Barrel Seal
 - the depressing of the plunger minimises back flow of blood, reducing the risk of blood exposure
- Plunger Head
 - when the plunger head is in line with the marker on the barrel, the guidewire is in line with the needle tip
- Large Flashback Chamber
 - the large transparent flashback chamber enables confirmation of vessel puncture
- Round Tip Guidewire
 - the round tipped guidewire aids smooth securing of the blood vessel
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- manufactured in BD Vialon™, a unique polyurethane, radiopaque biocompatible catheter material (please see page 48 for further information)



BD Instaflash™ Technology



Gauge Size	Outer Diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Insyte-A™ Arterial Catheter					
22 G	0,9	38	blue	20 / 80	382802
20 G	1,1	38	pink	20 / 80	382805

BD

- Catheter designed for insertion into peripheral arteries, for the determination of blood gas estimations or monitoring blood pressure
- 'Venflon-like' insertion characteristics
- Wings to facilitate fixation
- Class IIa according to European Medical Devices Directive 93/42/EEC
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked
- The Floswitch On/Off device prevents back-flow thereby reducing the risk of air emboli and potential contamination of hospital staff by exposure to blood



Gauge Size	Outer Diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Arterial Cannula					
20 G	1,1	45	red	25 / 500	682245

Arterial Catheters

BD Arterial Needle for Seldinger Technique

- Two-part arterial needle, specifically designed for the Seldinger Technique
 - Solid Inner Needle
 - precisely ground, extremely sharp, short bevel arterial point assures percutaneous entry with minimum trauma and vessel roll
 - solid design eliminates coring
 - Outer Thin-Wall Needle
 - close fitting taper provides smooth entry with minimum trauma
 - tapered, square-cut distal end of outer cannula facilitates smooth passage of guide wire through needle
 - lightweight polypropylene
- winged handle provides positive control
 - manufactured to exacting tolerances to accept guide wire with minimum blood leakage and associated clotting
 - funnel-shaped interior of outer needle hub facilitates insertion of guide wire
- Sterilization method is Ethylene Oxide
 - Sterile, single use
 - CE marked



Gauge Size	Outer Diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD™ Arterial Needle for Seldinger Technique					
18 G	1,2	73	pink	10 / 50	408273

BD M/P (Multi Purpose) Procedure Needle

- A multi purpose needle for traditional or digital angiography. Winged design provides positive control for arterial or venous punctures. Funnel shaped hub facilitates insertion of guide wire.
- Lightweight plastic BD Luer-Lok™ hub provides excellent "feel" and balance.
- Thin-wall needle for maximum I.D., minimum O.D. Sharp bevel for easy vessel entry
- single wall puncture (M/P) procedure needle
 - one part needle with wings
 - .032 in. - .038 in. guide wire
- Sterilization method is Ethylene Oxide
 - Sterile, single use
 - CE marked



Gauge Size	Outer Diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD M/P Procedure Needle					
18 G	1,2	70	pink	10 / 50	408295

Anaesthesia Systems

Spinal Anaesthesia

BD Whitacre Pencil Point Spinal Needle

BD Whitacre Needles are manufactured with proprietary needle tipping, side-port and unique low-glare hub technologies.

- Tip designed to spread the dural fibres and help reduce the occurrence of post dural puncture headache¹
- The thin wall design is intended for improved visualization of CSF
- Distinct “pop” as the pencil point penetrates the dura
- Pencil point design intended to provide better “tip strength” to minimize bending or breakage
- Precision-formed side hole allows / facilitates directional flow of anaesthetic agents
- Translucent hub for ease of CSF visualization

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked

¹ Vallejo MC, Mandell GL, Sabo DP, Ramanathan S (2000): Postdural puncture headache: a randomized comparison of five spinal needles in obstetric patients. *Anesth Analg.* 91(4):916-20



Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Whitacre Pencil Point Spinal Needle						
27 G	3 1/2"	0,4	90	grey	25 / 200	402051
27 G	4 11/16"	0,4	119	grey	10 / 50	409443
25 G	3 1/2"	0,5	90	blue	25 / 200	402050
25 G	4 11/16"	0,5	119	blue	10 / 50	409442
24 G	3 1/2"	0,55	90	purple	25 / 200	405104
22 G	3 1/2"	0,7	90	black	25 / 200	401995
BD Whitacre Pencil Point Spinal Needle – with integral introducer 22 G 1 1/4"						
27 G	3 1/2"	0,4	90	grey	25 / 200	405075
27 G	4 1/16"	0,4	103	grey	25 / 200	405113
BD Whitacre Pencil Point Spinal Needle – with integral introducer 20 G 1 1/4"						
25 G	3 1/2"	0,5	90	blue	25 / 200	405076
25 G	4 1/16"	0,5	103	blue	25 / 200	405112
Introducer						
20 G	1 1/4"	for 24/25/26 G Needle	32	transparent	25 / 200	405260
22 G	1 1/4"	for 27 G Needle	32	transparent	25 / 200	405261

Anaesthesia Systems

Spinal Anaesthesia

BD Spinal Needle Quincke Type Point

- Key/slot arrangement of stylet and cannula hubs facilitates aligned needle and stylet bevel orientation
- Translucent window hub (featuring contact clarity) that helps visualization of CSF
- Fitted stylet designed to reduce tissue coring
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Spinal Needles Quincke Type Point						
27 G	3 1/2"	0,4	90	grey	25 / 200	405259
26 G	3 1/2"	0,4	90	brown	25 / 200	405258
25 G	3 1/2"	0,5	90	orange	25 / 200	405257
25 G	3"	0,5	75	orange	25 / 200	405246
23 G	3 1/2"	0,6	90	turquoise	25 / 200	405240
22 G	3 1/2"	0,7	90	black	25 / 200	405256
22 G	3"	0,7	75	black	25 / 200	405255
22 G	2 1/2"	0,7	63	black	25 / 200	405244
22 G	1 1/2"	0,7	38	black	25 / 200	405254
20 G	3 1/2"	0,9	90	yellow	25 / 200	405253
20 G	3"	0,9	75	yellow	25 / 200	405252
20 G	1 1/2"	0,9	38	yellow	25 / 200	405251
19 G	3 1/2"	1,1	90	tan	25 / 200	405250
19 G	3"	1,1	75	tan	25 / 200	405249
18 G	3 1/2"	1,2	90	pink	25 / 200	405248
18 G	3"	1,2	75	pink	25 / 200	405247

Anaesthesia Systems

Spinal Anaesthesia

BD Spinal Needle Quincke Type Point

- Key/slot arrangement of stylet and cannula hubs facilitates aligned needle and stylet bevel orientation
- Translucent window hub (featuring contact clarity) that helps visualization of CSF
- Fitted stylet designed to reduce tissue coring
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Spinal Needle Quincke Type Point						
25 G	4 11/16"	0,5	119	blue	10 / 50	405234
22 G	7"	0,7	178	black	10 / 50	405149
22 G	5"	0,7	127	black	10 / 50	405148
20 G	6"	0,9	152	yellow	10 / 50	405211
18 G (intradiscal)	6"	1,2	152	pink	10 / 50	408360
BD Spinal Needle Quincke Type Point - with introducer 22 G 1 1/4"						
Spinal Needle 27 G	3 1/2"	0,4	90	grey	25 / 200	405069
BD Spinal Needle Quincke Type Point - with introducer 20 G 1 1/4"						
26 G	3 1/2"	0,45	90	brown	25 / 200	405065
25 G	3 1/2"	0,5	90	orange	25 / 200	405084
Introducers						
20 G	1 1/4"	for 24/25/26 G Needles	32	transparent	25 / 200	405260
22 G	1 1/4"	for 27 G Needles	32	transparent	25 / 200	405261
BD Neonatal Needle						
25 G (short bevel)	1"	0,5	25	blue	25 / 200	405243
25 G (short bevel)	2"	0,5	51	blue	25 / 200	405245

Perisafe Plus!™ Epidural Procedure Kit

A selection of kits for epidural anaesthesia procedures with BD Weiss needles

- BD Perisafe Modified Tuohy needle point designed for greater penetration characteristics
- Shorter point length than standard epidural needles designed to minimize accidental dural puncture
- Key/slot arrangement of stylet and cannula hubs facilitates aligned needle and stylet bevel orientation
- BD Epilor™ Loss-of-Resistance (L.O.R.) plastic syringe designed for better detection of epidural space

Perisafe Plus!™ Epidural Procedure Kits:

- Weiss Epidural Needle
- Epidural Catheter 90 cm
- Catheter Connector
- 0,22 µm Flat Filter

Standard components included in all Mini-Perisafe Plus!™ Epidural Procedure Kits:

- Same as Perisafe Plus!™, however does not include L.O.R. Syringe
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Standard components included in all

Weiss Epidural Needle	L.O.R. Plastic Syringe	Open End Epidural Catheter	Closed End Epidural Catheter	Catheter Connector	0,22 µm Flat Filter	Shelf Box / Shipping Case	Reference Number
Perisafe Plus!™ Epidural Procedure Kit							
17 G	✓	–	19 G	✓	✓	12 / 12	401506M
18 G	✓	–	20 G	✓	✓	12 / 12	401508M
18 G	✓	20 G	–	✓	✓	12 / 12	401509M
Mini-Perisafe Plus!™ Epidural Procedure Kit							
18 G	–	–	20 G	✓	✓	12 / 12	401511M

Anaesthesia Systems

Epidural Products

Single Components for Epidural Anaesthesia

BD Epilor™

Loss-of-Resistance Syringe

BD Epilor™ Loss-of-Resistance Syringe (LOR) has a double ribbed stopper to prevent leakage during aspiration.

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



BD Tuohy and

BD Weiss Epidural Needles

- Epidural needles with locking base and 'modified Tuohy' bevel intended to minimize accidental dural puncture
- Needle graduation marks every 10 mm. Paediatric epidural needle graduation marks every 5 mm.
- Needle equipped with moulded plastic hub
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Epilor™ Loss-of-Resistance syringe						
7 ml L.O.R. Plastic Syringe with BD Luer-Lok™					10 / 50	405291
7 ml L.O.R. Plastic Syringe Luer Slip					10 / 50	405292
BD Tuohy Epidural needle without wings						
17 G	3 1/2"	1,5	90	grey	10 / 50	405092
18 G	3 1/2"	1,2	90	pink	10 / 50	405018
20 G	3 1/2"	0,9	90	yellow	10 / 50	405028
BD Weiss Epidural needle with wings						
16 G	3 1/2"	1,5	90	transparent	10 / 50	405194
17 G	3 1/2"	1,5	90	grey	10 / 50	405191
17 G	5"	1,5	127	grey	10 / 50	405193
18 G	3 1/2"	1,27	90	pink	10 / 50	405188
18 G	5"	1,27	127	pink	10 / 50	405190
20 G	3 1/2"	0,9	90	yellow	10 / 50	405186

Anaesthesia Systems

Combined Spinal / Epidural Anaesthesia Sets

BD Durasafe Plus™! with CSE lock

The BD range of Combined Spinal / Epidural (CSE) anaesthesia sets has been designed to help ensure correct placement of the epidural needle.

BD Durasafe Plus™! CSE set intended for rapid onset of anaesthesia, extended duration, and better pain management.

- CSE set system with Withacre needle length matched to epidural needle: spinal needle can extend up to 15 mm beyond epidural needle point
- Secure locking device permits the spinal needle to lock at any extension (0-15mm) distal to the epidural needle tip
 - Whitacre spinal needle can rotate 360° when the mechanism is locked
 - allows spinal needle to move freely through mechanism when unlocked
 - the locking device is designed for ease of use and stability for better user control
- Standard components included in all kits:
 - Weiss Epidural Needle
 - Loss-of-Resistance Syringe (L.O.R)
 - BD Luer-Lok™ cap
 - Catheter Connector
 - Round Tip Epidural Catheter with Introducer
 - Tuohy Borst Adaptor (not shown)
 - 0,22 µm Flat Filter
 - BD Whitacre Spinal Needle
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Weiss Epidural Needle	Spinal Needle	L.O.R. Plastic Syringe	Epidural Catheter	Shelf Box / Shipping Case	Reference Number
BD Durasafe Plus™! with CSE lock					
17 G	BD Whitacre 27 G	✓	19 G	6 / 30	400286
18 G	BD Whitacre 27 G	✓	20 G	6 / 30	400287

1

2

3

4

6

7

8

9

10

11

12


13

Surgical Accessories

Biopsy Needles

BD Chiba Biopsy Needle

- 22G (7/10 mm) thin wall, short bevel (26°) needle with accurate grip which reduces the risk of vessel and tissue damage. Translucent base aids visualisation of vascular breaks
- Hub locking system to prevent biopsy hazards
- Two different needle lengths available
- Needle graduated in centimetres to measure insertion depth
- Transhepatic cholangiography
- Chemonucleolysis:
 - the intradiscal needle is used as a guide for the Chiba needle the distal end of which has a 2 cm curve to improve intradiscal insertion
- Anaesthesiology:
 - Sciatic block
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Indications

- Soft tissue biopsy

Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Chiba Biopsy Needle						
22 G	8"	0,72	203	black	10 / 50	408269
22 G	6"	0,72	152	black	10 / 50	408268

BD Westcott Biopsy Needle

- For direct percutaneous fine needle aspiration biopsy
- Available in two gauges and three lengths for maximum versatility for both deep and superficial masses
- This needle has proven use for many biopsy applications including that of pulmonary, hilar and mediastinal masses, as well as lesions in the liver, pancreas and other areas
- The BD Westcott Biopsy needle is a 2-piece, thin-wall needle modified to include a slotted opening measuring 2.2 mm in length located approximately 3 mm from the needle tip
- The slot creates a second cutting edge, in addition to the needlepoint, enabling aspirations of larger amounts of material without any apparent increase in complications
- For histological, cytological and bacteriological sampling
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Gauge Size	Length (Inch)	Outer diameter (mm)	Length (mm)	Color Code	Shelf Box / Shipping Case	Reference Number
BD Westcott Biopsy needle						
22 G	7"	0,7	178	black	10 / 50	408267
22 G	5 1/2"	0,7	135	black	10 / 50	408266
22 G	3 1/2"	0,7	90	black	10 / 50	408265
20 G	7"	0,9	178	yellow	10 / 50	408264
20 G	5 1/2"	0,9	135	yellow	10 / 50	408263
20 G	3 1/2"	0,9	90	yellow	10 / 50	408262

Surgical Accessories

BD Heimlich Chest Drain, Brush/Sponge

BD Heimlich Chest Drain

- The Heimlich chest drain valve is a specially-designed flutter valve used to replace underwater bottles in chest drainage
- Constructed of rubber tubing, the valve is compressed at one end to form leaflets that control unidirectional flow
- The valve is encased in a transparent plastic chamber with tapered ends
- The proximal end of the casing can be attached to most chest catheters and the distal end to tubing that empties into a plastic bag
- If desired, the distal end can be connected to regulated suction
- PVC free, contains latex
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Shelf Box / Shipping Case	Reference Number
BD Heimlich Chest Drain	10 / 10	373460

BD E-Z Scrub™ Brush/Sponge

BD single use scrub brush has soft polyethylene bristles on one side and a polyether sponge on the other. Ergonomically designed for greater user comfort and effectiveness.


Rigidity

- The standard rigidity of BD E-Z Scrub™ scrub brush has been designed to provide the efficiency and comfort needed for surgical scrub brushes. The tapered design allows for two degrees of softness, appropriate for different areas of skin conditions. The inner recessed bristles are stiffer and designed specially

to remove dirt from beneath the nails. The sponge provides cleaning for gentle use on irritated areas, or when/where severe mechanical action or brushing is not required, such as arms and forearms.

Recommended use

- Surgical scrubbing of hands before wearing sterile gloves for aseptic procedures
- Sterilization method is Radiation
- Sterile, single use



Beschreibung	Shelf Box / Shipping Case	Reference Number
BD E-Z Scrub™ brush / sponge		
without active ingredient, dry, sterile	30 / 300	371607

1

2

3

4

5

7

8

9

10

11

12

13

Sharps Disposal Systems

BD™ Sharps Disposal Systems

BD™ Sharps Disposal Systems - focused on "Safety".

The new range of BD™ Sharps Disposal Containers is designed to optimise "Point of Use" disposal to help reduce the risks of blood exposure and needle stick injuries associated with the disposal process.

BD™ Sharps Container come in sizes of 1,5 l - 3,0 l - 5,0 l - 7,0 l

- Functional ease of use - visible and audible locking system to ensure proper assembly
- Visible temporary and permanent closure features
- Temporary and permanent closure engagement without touching the access area
- Retention valve against spillage due to optimal container opening for increased ease of disposal while restricting access for safety
- Manually inviolable when container is permanently locked
- Clear visualisation to the maximum fill line
- Rounded corners to ensure optimal impact resistance if dropped to avoid cracking
- Accompanied by a range of brackets and accessories focusing on "correct placements" and "secure fitment"¹
- BD™ Sharps Containers are made of polypropylene and may be incinerated and autoclaved

Regulatory Compliance

- Compliant with the British Standard BS 7320 : 1990²
- Compliant with the French Standard NFX 30 - 500³
- Compliant with the UN Standard, UN 3291⁴

¹ "Selecting, Evaluating and Using Sharps Disposal Containers" NIOSH/CDC Publication No. 97-111

² not for BD™ Sharps Container 0.87 l

³ not for BD™ Sharps Containers 0.45 and 22.7 l

⁴ not for BD™ Sharps Containers 0.87 and 22.7 l



Nominal Capacity (l)	Usable Capacity (l)	Dimensions (cm)	Shelf Box / Shipping Case	Reference Number
BD™ Sharps Disposal Container				
0,45	0,36	10,5 x 5,0 x 17,0	100 / 100	302434
0,87	0,7	8,9 x 7,4 x 16,8	50 / 50	367202
1,5	1,2	16,5 x 16,0 x 20,5	40 / 40	305624
3,0	2,4	18,2 x 18,1 x 26,0	25 / 25	305625
5,0	4,0	30,4 x 22,2 x 20,8	20 / 20	305626
7,0	5,6	30,4 x 22,2 x 25,5	16 / 16	305627
22,7	18,2	23,0 x 22,7 x 44,5	12 / 12	300479

Sharps Disposal Systems

Brackets and Accessories

Using brackets with the correctly selected sharps container is critical to ensure the correct placement and the secure fitment of the sharps container at the Point of Use for maximum safety.

Clearly identifying a place for the sharps container using a bracket will help to ensure the availability in the appropriate setting at all times.

The secure fitment of the sharps container in a bracket ensures maximum stability and allows the attention to be focused on the patient and not on the disposal process.

Desktop Brackets (fitting for 1,5 - 7 l BD™ Sharps Container)

- For placement on flat surfaces or secured to walls
- Attachable to rails with the addition of a Universal Rail Adapter
 - Universal Rail Adapter comes in 2 sizes adjustable to round or square rails

Mounting accessories for all Desktop Bracket sizes

- on a smooth horizontal surface
- on a vertical wall
- on round and square rails
- on a medication cart

Hanging Brackets (fitting for 1,5 - 7 l BD™ Sharps Container)

- for a placement on medication carts

The BD Desktop and Hanging Brackets are made of metal with a plastic coating to ensure that they are "easy clean" and remain "rust free" for safe hygiene practices and optimum use conditions.



Description	Secure fitment of							Shelf Box / Shipping Case	Reference Number
	0,45 l	0,87 l	1,5 l	3,0 l	5,0 l	7,0 l	22,7 l		
Brackets supports and accessories									
Plastic bracket for placement on horizontal surfaces	✓	–	–	–	–	–	–	12 / 12	305964
Metal cup bracket for placement on horizontal surfaces	–	✓	–	–	–	–	–	50 / 50	305970
Desktop bracket	–	–	✓	–	–	–	–	6 / 6	305645
Desktop bracket	–	–	–	✓	–	–	–	6 / 6	305250
Desktop bracket	–	–	–	–	✓	–	–	6 / 6	305253
Desktop bracket	–	–	–	–	–	✓	–	6 / 6	305255
Hanging bracket	–	–	✓	–	–	–	–	6 / 6	305710
Hanging bracket	–	–	–	✓	–	–	–	6 / 6	305251
Hanging bracket	–	–	–	–	✓	–	–	6 / 6	305252
Hanging bracket	–	–	–	–	–	✓	–	6 / 6	305254
Accessories									
Universal Adapter (small)	fitting for all Desktop brackets							6 / 6	305636
Universal Adapter (large)	fitting for all Desktop brackets							6 / 6	305637
Replacement Suction Cups	fitting for all Desktop brackets							12 / 144	305968

1

2

3

4

5

6

7

8

9

10

11

12

13

Closed System Drug Transfer Device

Take a closer look at the first FDA cleared CSTD indicated to reduce exposure to hazardous drugs

BD PhaSeal™ Closed System Drug Transfer Device is the first CSTD cleared under the FDA's newly created product code, ONB: Closed Antineoplastic and Hazardous Drug Reconstitution and Transfer System. Based on:

- Nothing Out –
No escape of hazardous drug or vapor concentration (leakproof & airtight)
- Nothing In –
No transfer of environmental contaminants
- Prevention of Microbial Ingress

New FDA code clarifies misconceptions in the market place*

- Establishes a specific product code for CSTD's indicated to reduce exposure to hazardous drugs
- Requires data to prove a system is closed for use with hazardous drugs
- Requires data to demonstrate reduction in healthcare worker exposure

New

ONB: Closed Antineoplastic and Hazardous Drug Reconstitution and Transfer System

VS

Existing

LHI: Set, I.V. Fluid Transfer
FPA: Intravascular Administration Set
FMF: Piston Syringe

	Specific to CSTDs	Requires Proof the system is closed	Requires Data to prove HCW exposure reduction
ONB	Yes	Yes	Yes
LHI / FPA / FMF	No	No [†]	No

* Data on file with BD

† Unless system is expressly indicated as a closed system drug transfer device (CSTD) that mechanically prohibits the transfer of environmental contaminants into the system and the escape of drug or vapor concentrations outside the system.

Closed System Drug Transfer Device

Protector

BD PhaSeal™ Protector – Drug Vial Access Device

The Protector is a pressure equalization device that is permanently attached to the vial and is used in the preparation of either powder or liquid drugs.

It provides airtight safety due to the sealed expansion chamber and dry leak-proof connections.

The Protector establishes and maintains neutral pressure when air or fluid is injected into or aspirated from the vial.

- For drug vials with neck diameters from 13 to 32 mm

- Reference Number 515103 and 515106 are packed in multipacks of 5 units each

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	For drug vials with neck diameter of (mm)	Pressure Equalization (ml)	Color Code	Shelf Box / Shipping Case	Reference Number
BD PhaSeal™ Protector					
Protector (P14)	13	20	blue	50 / 200	515100
Protector (P21)	20	20	red	50 / 200	515102
Protector (P21 Multi)	20	20	red	250 case	515103
Protector (P28)	28	60	white	25 / 100	515104
Protector (P50)	20	50	green	25 / 100	515105
Protector (P50 Multi)	20	50	green	150 case	515106
Protector (P53)	32	50	purple	25 / 100	515107

1

2

3

4

5

6

7

8

9

10

11

12

13

Closed System Drug Transfer Device

Injector

BD PhaSeal™ Injector – Syringe Safety Device

The Injector Luer Lock attaches to a standard syringe or IV line for a dry, leakproof connection during drug preparation and administration. It ensures a closed drug transfer by means of tightly sealed, double elastomeric membranes.

- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Injector (N30C)

- Injector (N30C) transfer device that attaches to a disposable syringe with standard Luer connection
 - with red protection cap

Injector (N35 and N35C)

- Injector (N35) and Injector (N35C) transfer devices that attach to a disposable syringe or IV tubing with a standard Luer Lock connection
 - Injector (N35) without protection cap
 - Injector (N35C) with clear protection cap
- Reference Number 515005 is packed in multipacks of 5 units each

Description	Protection cap	Shelf Box / Shipping Case	Reference Number.
BD PhaSeal™ Injector			
Injector Luer (N30C) for syringes with Luer connection	red	50 / 200	515001
Injector Luer Lock (N35) for syringes or IV bag with Luer Lock connection	–	50 / 200	515003
Injector Luer Lock (N35C) for syringes or IV bag with Luer Lock connection	clear	50 / 200	515004
Injector Luer Lock (N35C Multi) for syringes or IV bag with Luer Lock connection	clear	250 case	515005

Closed System Drug Transfer Device

Connector

BD PhaSeal™ Connector – For Closed Administration of Hazardous Drugs

The Injector Luer Lock attaches to a standard syringe or IV line for a dry, leakproof connection during drug preparation and administration.

It ensures a closed drug transfer by means of tightly sealed, double elastomeric membranes.

Connector Luer Lock (C35)

- Transfer device that attaches to a disposable syringe with standard Luer connection, with cap
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Connector Luer Lock (C45)

- Attaches to the patient's IV line, providing a closed connection to the Injector
- Standard Luer Lock connection
- This Connector is designed for use with needle-free ports with internal pistons
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Y-Site Connector (C80)

- Extension tubing with built-in BD PhaSeal™ Connector
- Total length: 21 cm
- Priming volume: 1,6 ml
- Contains PVC (DEHP free)
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Length (cm)	Priming Volume (ml)	Shelf Box / Shipping Case	Reference Number
BD PhaSeal™ Connector				
Connector Luer Lock (C35)			50 / 200	515200
Connector Luer Lock for use with needle-free ports with internal pistons (C45)			50 / 200	515202
Y-site connector extension tubing with built-in BD PhaSeal™ Connector (C80)	22	1,7	30 / 120	515304

Closed System Drug Transfer Device

Administration products – IV bag access

BD PhaSeal™ Administration products – IV bag access

Infusion Set (C50)

- Standard IV infusion set with built-in PhaSeal Connector for closed transfer to IV bag
- Total length: 188 cm
- Priming volume: 17,5 ml
- Inline filter: 15 µm
- Drops/ml: 20
- Contains PVC (DEHP free)
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Secondary Set (C60), Drip Chamber

- Secondary set with built-in PhaSeal Connector for closed transfer to IV bag.
- Total length: 102 cm
- Priming volume: 12,5 ml
- Drops/ml: 20
- Contains PVC (DEHP free)
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Secondary Set (C61)

- Secondary set with in-built PhaSeal Connector for closed transfer to IV bag
- Total length: 49 cm
- Priming volume: 3 ml
- Contains PVC (DEHP free)
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Length (cm)	Priming Volume (ml)	Inline Filter	Drops (ml)	Shelf Box / Shipping Case	Reference Number
BD PhaSeal™ administration products - IV bag access						
Infusion Set (C50)	188	17	15 µm	20	80 / 80	515300
Secondary Set Drip Chamber (C60)	102	12,5	–	20	20 / 80	515301
Secondary Set (C61)	49	3,0	–		30 / 120	515302

Closed System Drug Transfer Device

Administration products – IV bag access

BD PhaSeal Administration products – IV bag access

Infusion Adapter (C70)

- Connection between IV bag and non-vented infusion set
- Built-in PhaSeal Connector for closed addition to IV bag
- Contains PVC and DEHP
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



L Connector (C90)

- The L Connector provides a dry connection to semi-rigid plastic containers or glass bottles for closed addition of hazardous drugs
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Infusion Adapter (C100)

- connection between IV bag and non-vented infusion set. Built-in BD PhaSeal™ Connector for closed transfer to IV bag
- Reference Number 515307 is packed in multipacks of 5 units each
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Cyto Admin 4 Line Set (C110)

- The benefits of the female Luer-Lock adapters with integrated check valves also allow applications where medication is to be administered by an infusion pump
- Contains PVC (DEHP free) and TOTM plasticizer
- The vented spike with check valve also permits use of bottles or bags
- Sterilization method is Ethylene Oxide
- Sterile, single use
- CE marked



Description	Length (cm)	Shelf Box / Shipping Case	Reference Number
BD PhaSeal™ administration products - IV bag access			
Infusion Adapter (C70)	17	50 / 50	515303
L Connector (C90)		50 / 200	515305
Infusion Adapter (C100)		50 / 200	515306
Infusion Adapter (C100), Multipack		250 case	515307
Cyto-Ad Z®*		60	515308

* Registered trademark of Codan

Closed System Drug Transfer Device

Application Products

BD PhaSeal™ Application products


Infusion Clamp (M25)

- Clamp for securing the connection of Injector Luer Lock (N35) to the Connector
- Non-sterile
- CE marked



Assembly Fixture (M12)

- Mechanical device used for consistent and simple attachment of a Protector to a drug vial
- Non-sterile
- CE marked



Description	Shelf Box / Shipping Case	Reference Number
BD PhaSeal™ Application Products		
Infusion Clamp (M25)	25 / 100	515403
Assembly Fixture (M12)	1 / 1	515400

1

2

3

4

5

6

7

8

9

10

11

12

13

Social Responsibility - Trusted Partners

Social Responsibility - Trusted Partners

The BD “Trusted Partners” campaign was launched in early 2003. The campaign highlights the partnerships BD has formed with philanthropic, government and nongovernmental agencies for more than 100 years in order to make a significant impact on healthcare and to save lives throughout the world. The “Trusted Partners” initiative focuses on the relationship of BD with prominent nonprofit organizations and describes the work being done through the partnerships to address some of the world’s most important and fundamental health issues. Through philanthropy, sustainable business efforts and the development of pioneering technologies, BD is helping to support the development of an effective AIDS vaccine, eliminate deadly diseases (including measles and maternal and neonatal tetanus), protect against smallpox, fight the spread of diseases and provide relief in the wake of natural and manmade disasters. These initiatives are an essential element of the commitment of BD to pursue its purpose of “Helping all people live healthy lives.”



DIABETES:

The leader in diabetes care since 1924, BD has the world’s most complete array of insulin delivery and monitoring products, donates millions of dollars to research and develops educational programs for diabetes patients.



AMERICARES:

Since 1986, BD and AmeriCares have partnered to make a positive impact on healthcare globally by responding to emergency medical needs and supporting humanitarian assistance programs.



HIV/AIDS:

BD is committed to the battle against HIV/AIDS for future generations on many fronts with partner organizations, such as the William J. Clinton Presidential Foundation, the Global Business Coalition on HIV/AIDS, the International AIDS Vaccine Initiative and Save the Children.



HEALTHCARE WORKER SAFETY:

In 1991, BD was founding sponsor of the University of Virginia International Healthcare Worker Safety Center and primary supporter of the EPINet™ Surveillance system, a database that tracks occupational blood exposures. This partnership demonstrates BD’s leading efforts to protect caregivers by designing devices with safety features.

Social Responsibility - Trusted Partners

Social Responsibility - Trusted Partners



IAVI:

BD has partnered with the International AIDS Vaccine Initiative (IAVI) to evaluate and identify vaccines that can prevent development of the AIDS virus.



PROJECT HOPE:

Through an expansive "Train the Trainer" program, BD is performing an influential role in the education of Chinese people about the prevention and treatment of diabetes.



UNICEF:

The purpose of the UNICEF program is to eliminate maternal and neonatal tetanus - a deadly health threat in many developing nations - by the year 2005.



AMERICAN RED CROSS:

BD's partnership with the American Red Cross supports the Measles Initiative, a commitment to vaccinate 200 million African children over 5 years.



BIOTERROR:

BD was on the frontline for polio vaccines and once again, BD is at the forefront in protection against smallpox.

1

2

3

4

5

6

7

8

9

10

11

12

13

Annex

General Conditions of Sales

BECTON DICKINSON INTERNATIONAL
Branch of BECTON DICKINSON B.V. ("BDI")
Erembodegem Dorp 86
9320 - Erembodegem (Belgium)

1. General

- 1.1 The General Sales Conditions shall apply to all offers (including offers by means of pro forma invoices), acceptances of purchase orders, sales, and deliveries of BDI.
- 1.2 The General Sales Conditions shall be deemed accepted by the Buyer if Buyer has not in writing clearly objected to them.
- 1.3 BDI objects to any general terms and conditions of the Buyer.
- 1.4 All sales shall be non-exclusive and shall not confer any distribution or agency rights unless otherwise specifically agreed.

2. Offers, oral statements and amendments

- 2.1 All offers by BDI shall be subject to alterations without prior notice. Any offer of BDI shall be deemed withdrawn if it has not been accepted within the period there mentioned or within fourteen (14) days after its issue.
- 2.2 Any oral statement, addition or amendment to the General Sales Conditions and/or the documents of the related transaction shall only become valid if confirmed by BDI in writing.

3. Purchase orders

- 3.1 Any purchase order shall be subject to BDI's approval but, however, may be accepted impliedly by simply effecting delivery.
- 3.2 Any purchase order shall identify the country of final destination.
- 3.3 Products may normally be expected to be available, ready for shipment, within reasonable time from BDI's acceptance of a purchase order but delivery, however, remain subject to BDI's manufacturing capabilities.
- 3.4 Any purchase order shall be of a minimum amount of 2,500.- US dollars. The Buyer placing a purchase order below such minimum value shall pay to BDI a handling charge of 250.- US dollars.

4. Prices

- 4.1 All quotations of prices are in US \$ except as otherwise provided for in BDI's offer plus, even if not so mentioned in the respective offer, taxes or other charges as provided for by the applicable law.
- 4.2 The price payable by the Buyer for the products delivered shall be BDI's current price list which is effective at the date of delivery.

5. Payment

- 5.1 Payment shall become due at the date shown on the invoice or, if not mentioned, within 30 days net of the invoice's date.

5.2 Delayed payments shall bear interest of 1,5 % per month which shall be paid by the Buyer to BDI and Buyer shall additionally reimburse BDI all costs and damages associated with delayed payment, legal expenses included.

5.3 In no event shall the Buyer have the right to retain or set off payments unless BDI consented thereto in writing or a final and binding court decision confirms Buyer's claim.

5.4 If payment shall be effected by a letter of credit, BDI's Letter of Credit Policy shall apply.

6. Delivery and risk

6.1 Any delivery of BDI shall be "ex-works" (Incoterms 2010). However, on Buyer's request, BDI may provide, at Buyer's exclusive risk and cost, for transportation, legalization, inspection, and any handling charges on usual terms.

6.2 If BDI has notified Buyer that the Products ordered are ready to ship, Buyer shall take delivery without delay failing which all cost for storage and risk of loss or damage shall be for Buyer's account.

6.3 Products are normally available within the time frame given in the purchase order as accepted but, however, remain subject to manufacturing capabilities and schedules as well as transportation facilities.

6.4 On BDI's request Buyer shall furnish customs clearance or any other documentation pertinent to the transaction.

7. Warranty

7.1 All products delivered shall be in compliance with the applicable specifications for the respective period of their shelf life.

7.2 In case of products being instruments the warranty shall be twelve (12) months from the installation date at the ultimate customer's facility but in no event to exceed eighteen (18) months as of delivery to the Buyer. However, any instrument warranty shall only be given on the condition that the respective warranty and installation certificate of Becton was properly completed by Buyer.

7.3 The warranty shall neither cover products damaged in transit nor when the Buyer is in default of acceptance. Mishandling, repairs or alterations of products without BDI's authorization shall render any warranty void.

7.4 Any defects of products which would be apparent on reasonable inspection must be notified to BDI in writing within 72 (seventy-two) hours of reception of them by the Buyer failing which the products shall irrevocably be deemed accepted by the Buyer.



















7.5 Save as expressly provided for in this clause BDI shall have no liability or responsibility to the Buyer in respect of defective products and all warranties regarding the quality, condition, durability, merchantability, fitness for purpose, conformity with samples, or likewise are hereby expressly excluded.

BECTON DICKINSON INTERNATIONAL
Branch of BECTON DICKINSON B.V. ("BDI")
Erembodegem Dorp 86
9320 - Erembodegem (Belgium)

- | | |
|---|--|
| <p>7.6 Any warranty claim shall be in writing, using BDI's applicable complaint form, and identify the defects and their presumable cause.</p> <p>8. Retention of title
BDI reserves title to all products supplied to Buyer until such time as payment for the relevant delivery has been fully received by BDI.</p> <p>9. Liability
9.1 Notwithstanding mandatory provisions of any applicable law, BDI shall neither be liable for indirect or consequential damages nor for any other damages unless the latter were caused by BDI's gross negligence or willfulness.
9.2 BDI shall in no event be liable for any damage if it cannot be excluded that a fault of the Buyer or of a person for whom the Buyer is responsible caused the damage including, not limited to events identified in section 7.3 above.</p> <p>10. Force Majeure
Neither party shall be liable for any breach of any agreement governed by the General Sales Conditions to the extent that this is caused by Force Majeure circumstances i.e. events which are beyond the reasonable control of the affected party, such as, but not limited to, acts of god, strikes, work to rule, governmental acts, fire, flood, bad weather conditions, war, shortfall of raw materials, energy, transportation, regardless of whether these were unforeseeable or not.</p> <p>11. Termination
Any agreement governed by the General Sales Conditions may be terminated by BDI with immediate effect if:</p> <ul style="list-style-type: none"> (i) Without requirement for BDI to any further notice, the Buyer fails to pay at due dates, or (ii) the Buyer is in material breach of any other term or condition of the General Sales Conditions including the related agreement and has not remedied such breach within fifteen (15) days of having been notified of such breach, or (iii) The Buyer shall be or become bankrupt or insolvent or if there are instituted by or against it proceedings in bankruptcy or under insolvency laws or for its reorganization, receivership, liquidation, dissolution or change of ownership. <p>12. Amendments
The General Sales Conditions including any related agreement shall not be changed or modified except by written agreement duly signed by BDI and Buyer.</p> | <p>13. Applicable Law and Jurisdiction
The General Sales Conditions and any related transaction shall be governed and construed in accordance with the material law of Belgium except its choice of law principles and the United Nations Convention on Contracts for the International Sale of Goods of April 11, 1980 and the exclusive venue of any dispute relating to the General Sales Conditions or the transactions contemplated herein shall be the Courts of Aalst.</p> <p>14. Anti-boycott legislation
Any agreement referring to the General Sales Conditions in conflict with US anti-boycott or any other applicable anti-boycott legislation shall be resolved first by doing the necessary to conform with the applicable anti-boycott provisions, including, not limited to, amending such agreement and, if that does not resolve the conflict, the violating agreement shall be null and void as if it had never come into existence.</p> <p>15. Severability
Should any part, term or provision of the General Sales Conditions be declared by any court to be or be accepted by the parties as being in conflict with the law or unenforceable the validity and enforceability of the remainder of the General Sales Conditions shall not be affected thereby. The offending part or provision shall be replaced by an enforceable and valid part or provision of which the economic intent is as far as possible the same.</p> <p>16. Trade Compliance
BD and 'customer' agree to comply fully with all laws governing trans-border sales, re-sales, shipments and transfers of Products, Original Product or End Product, including but not limited to the export control laws of the United States and of any other countries, to the extent that such laws are applicable to the parties or the Product, Original Product or End Product. In this regard:</p> <ul style="list-style-type: none"> - BD's obligation to supply Product or Original Product to 'customer' is contingent upon receiving any required governmental authorizations. Upon BD's request, 'customer' agrees that it will furnish all information and documentation necessary for BD to obtain such authorizations. - 'Customer' agrees to comply with all conditions and destination control statements set forth on the invoice, bill of lading or other documents accompanying the shipment of BD Products or Original Products. - 'Customer' affirms that it will not knowingly use, resell or distribute Product, Original Product or End Product, directly or indirectly, for the development, production or proliferation of weapons of mass destruction (nuclear, chemical, or biological) or missile delivery systems, and/or for terrorist activities. |
|---|--|

2
3
4
5
6
7
8
9
10
11
12
13

Symbol Guidance

 REF Catalogue Number / Reference Number	 Use By	 Do Not Reuse
 CE Conformité Européenne	 Date of Manufacture	 Keep Dry
 LOT Batch Code	 Manufacturer	 Peel Here
 EC REP Authorized Representative in the European Community	 Do Not Use if Package is damaged	 Caution
 STERILE EO Sterilized Using Ethylene Oxide	 Keep Away from Sunlight	 Consult Instructions For Use
 STERILE R Sterilized Using Radiation	 STERILE Method of Sterilization Using Steam or Dry Heat	 Fragile, Handle with Care

BD Medical - Medical Surgical Systems products are designed and manufactured to meet applicable international Quality Management System requirements in accordance with the European Medical Devices Directive 93/42/EEC (MDD) and FDA QSR where applicable.

Hence, all BD manufacturing sites hold a valid certification against the applicable international Quality Management Systems standards, such as e.g. ISO 13485:2012, EN ISO 9001:2008.

Likewise, all BD manufacturing sites are subject of periodic external as well as internal audits to ensure compliance with aforementioned Quality Management System requirements.

Index

Numerical Index

300013.....	16	301355.....	16	305273.....	11
300014.....	16	301500.....	5	305274.....	11
300015.....	16	301504.....	14	305275.....	11
300026.....	16	301700.....	5	305283.....	11
300094.....	5	301750.....	5	305307.....	23
300136.....	17	301900.....	5	305309.....	23
300137.....	17	302187.....	15	305310.....	11
300138.....	17	302188.....	15	305311.....	11
300139.....	17	302200.....	5	305312.....	11
300185.....	15	302245.....	22	305313.....	11
300296.....	14	302247.....	22	305501.....	16
300300.....	5	302248.....	22	305502.....	16
300318.....	14	302249.....	22	305624.....	60
300320.....	14	302250.....	22	305625.....	60
300321.....	14	302251.....	22	305626.....	60
300322.....	14	302434.....	60	305627.....	60
300323.....	14	302436.....	8	305636.....	61
300330.....	14	302437.....	8	305637.....	61
300334.....	16	302834.....	29	305645.....	61
300400.....	5	303800.....	5	305710.....	61
300479.....	60	304000.....	5	305760.....	8
300506.....	14	304100.....	5	305770.....	8
300512.....	14	304300.....	5	305771.....	8
300513.....	14	304432.....	5	305781.....	8
300516.....	14	304434.....	5	305782.....	8
300600.....	5	304622.....	5	305783.....	8
300605.....	17	304727.....	5	305784.....	8
300613.....	17	304827.....	5	305785.....	8
300626.....	14	305180.....	7	305786.....	8
300629.....	17	305181.....	7	305833.....	11
300635.....	5	305207.....	18	305886.....	8
300637.....	5	305208.....	18	305887.....	8
300700.....	5	305209.....	18	305888.....	8
300800.....	5	305210.....	18	305889.....	8
300865.....	17	305211.....	7	305891.....	8
300866.....	17	305217.....	18	305892.....	8
300867.....	17	305218.....	18	305894.....	8
300869.....	17	305219.....	18	305895.....	8
300900.....	5	305220.....	18	305899.....	8
300912.....	15	305250.....	61	305900.....	9
300928.....	14	305251.....	61	305901.....	9
301000.....	5	305252.....	61	305902.....	9
301001.....	14	305253.....	61	305904.....	9
301155.....	5	305254.....	61	305905.....	9
301156.....	5	305255.....	61	305907.....	9
301190.....	17	305269.....	11	305915.....	9
301229.....	17	305270.....	11	305916.....	9
301231.....	17	305271.....	11	305917.....	9
301300.....	5	305272.....	11	305918.....	9

1

2

3

4

5

6

7

8

9

10

11

12

13

Index

Numerical Index

305921.....	9	381037.....	37	381944.....	36, 45
305945.....	10	381044.....	37	381947.....	36, 45
305946.....	10	381047.....	37	381954.....	36
305964.....	61	381054.....	37	381957.....	36, 45
305968.....	61	381057.....	37	382258.....	42
305970.....	61	381112.....	42	382259.....	42
306580.....	48	381123.....	42	382268.....	42
306581.....	48	381134.....	42	382269.....	42
306582.....	48	381137.....	42	382277.....	42
306583.....	48	381144.....	42	382287.....	42
306584.....	48	381147.....	42	382802.....	50
306585.....	48	381154.....	42	382805.....	50
307727.....	13	381157.....	42	382912.....	37
307728.....	13	381164.....	42	382923.....	37
307729.....	13	381167.....	42	382933.....	37
307730.....	13	381211.....	41	382934.....	37
307731.....	13	381212.....	41, 45	382937.....	37
307732.....	13	381223.....	41, 45	382944.....	37
307733.....	13	381234.....	41, 45	382947.....	37
307734.....	13	381237.....	41, 45	382954.....	37
307735.....	13	381244.....	41, 45	382957.....	37
307736.....	13	381247.....	41, 45	383318.....	38
307737.....	13	381257.....	41, 45	383319.....	38
307738.....	13	381267.....	41, 45	383328.....	38
307739.....	13	381311.....	41	383329.....	38
307740.....	13	381312.....	41, 45	383338.....	38
307741.....	13	381323.....	41, 45	383339.....	38
307742.....	13	381334.....	41, 45	383348.....	38
307743.....	13	381337.....	41, 45	383510.....	33
307960.....	21	381344.....	41, 45	383511.....	33
307961.....	21	381347.....	41, 45	383512.....	33
307962.....	21	381357.....	41, 45	383516.....	33
307964.....	21	381811.....	36	383517.....	33
307970.....	21	381812.....	36, 45	383518.....	33
307971.....	21	381823.....	36, 45	383519.....	33
307980.....	21	381833.....	36	383520.....	33
307981.....	21	381834.....	36, 45	383530.....	33
309050.....	14	381837.....	36, 45	383531.....	33
309110.....	14	381844.....	36, 45	383532.....	33
309628.....	16	381847.....	36, 45	383536.....	33
309649.....	15	381854.....	36	383537.....	33
309658.....	15	381857.....	36, 45	383539.....	33
367202.....	60	381867.....	36, 45	383540.....	33
371607.....	59	381911.....	36	383640.....	33
373460.....	59	381912.....	36, 45	383641.....	33
381012.....	37	381923.....	36, 45	383642.....	33
381023.....	37	381933.....	36	383646.....	33
381033.....	37	381934.....	36, 45	383647.....	33
381034.....	37	381937.....	36, 45	383648.....	33

Index

Numerical Index

383649.....	33	393224.....	35, 44	396368.....	25
383650.....	33	393226.....	35, 44	396370.....	25
383660.....	33	393227.....	35, 44	396371.....	25
383661.....	33	393228.....	35, 44	396373.....	25
383662.....	33	393229.....	35, 44	396377.....	25
383666.....	33	393230.....	35, 44	396390.....	25
383667.....	33	394074.....	29	396392.....	25
383668.....	33	394075.....	29	397015.....	29
383669.....	33	394080.....	29	397019.....	30
383670.....	33	394090.....	29	397251.....	24
383681.....	33	394251.....	44	397400.....	24
383682.....	33	394252.....	44	397509.....	24
383686.....	33	394253.....	44	397590.....	24
383687.....	33	394254.....	44	397608.....	24
383689.....	33	394255.....	44	397610.....	24
383691.....	34	394256.....	44	397611.....	24
383692.....	34	394257.....	44	397994.....	24
383693.....	34	394280.....	45	397997.....	24
383694.....	34	394281.....	45	400286.....	57
383695.....	34	394282.....	45	400287.....	57
385100.....	26	394283.....	45	401506M.....	55
385101.....	26	394284.....	45	401508M.....	55
385102.....	26	394285.....	45	401509M.....	55
385108.....	27	394286.....	45	401511M.....	55
385161.....	26	394287.....	45	401995.....	52
385162.....	26	394501.....	27	402050.....	52
385163.....	26	394600.....	28	402051.....	52
385164.....	26	394601.....	28	405018.....	56
386411.....	31	394602.....	28	405028.....	56
386412.....	31	394605.....	28	405065.....	54
391349.....	43	394926.....	28	405069.....	54
391350.....	43	394936.....	28	405075.....	52
391451.....	40	394945.....	28	405076.....	52
391452.....	40	394951.....	28	405084.....	54
391453.....	40, 44	394961.....	28	405092.....	56
391454.....	40, 44	394971.....	28	405104.....	52
391455.....	40, 44	394982.....	28	405112.....	52
391456.....	40, 44	394983.....	28	405113.....	52
391457.....	40, 44	394995.....	28	405148.....	54
392015.....	43	394997.....	28	405149.....	54
392020.....	43	396350.....	25	405186.....	56
393202.....	39, 44	396351.....	25	405188.....	56
393204.....	39, 44	396352.....	25	405190.....	56
393206.....	39, 44	396353.....	25	405191.....	56
393207.....	39, 44	396354.....	25	405193.....	56
393208.....	39, 44	396357.....	25	405194.....	56
393209.....	39, 44	396359.....	25	405211.....	54
393210.....	39, 44	396361.....	25	405234.....	54
393222.....	35, 44	396363.....	25	405240.....	53

1

2

3

4

5

6

7

8

9

10

11

12

13

Index

Numerical Index

405243.....	54	515302.....	66
405244.....	53	515303.....	67
405245.....	54	515304.....	65
405246.....	53	515305.....	67
405247.....	53	515306.....	67
405248.....	53	515307.....	67
405249.....	53	515308.....	67
405250.....	53	515403.....	68
405251.....	53	515404.....	68
405252.....	53	682245.....	50
405253.....	53		
405254.....	53		
405255.....	53		
405256.....	53		
405257.....	53		
405258.....	53		
405259.....	53		
405260.....	54		
405261.....	54		
405291.....	56		
405292.....	56		
408262.....	58		
408263.....	58		
408264.....	58		
408265.....	58		
408266.....	58		
408267.....	58		
408268.....	58		
408269.....	58		
408273.....	51		
408295.....	51		
408360.....	54		
409442.....	52		
409443.....	52		
515001.....	64		
515003.....	64		
515004.....	64		
515005.....	64		
515100.....	63		
515102.....	63		
515103.....	63		
515104.....	63		
515105.....	63		
515106.....	63		
515107.....	63		
515200.....	65		
515202.....	65		
515300.....	66		
515301.....	66		

Index

Alphabetical Index

A

Acacia Inc. Specialty Extension Sets with BD Q-Syte™	26
Adapter	
BD Q-Syte™ Vial Access Adapter	27
Infusion Adapter, BD PhaSeal™	67
Anaesthesia Systems	52–57
Combined Spinal / Epidural Anaesthesia Sets	57
Epidural Products	55
Spinal Anaesthesia	52
Angiocath™	42
Arterial Cannula	50
Arterial Catheters	50–51
Arterial Needle for Seldinger Technique	51
Autodisable Immunization Syringes	22
BD SoloShot™ Mini	22

B

BD Cannulae / Catheter Conversions	47
BD PhaSeal™	
Connector	65
Injector	64
Protector	63
Biopsy Needles	
BD™ Westcott Biopsy Needle	58
Blunt Fill Needle	7

C

Caps and Plugs	29
Chiba Biopsy Needle	58
Closed IV Catheter System	
BD Nexiva™	32, 33
BD Nexiva™ Diffusics™	34

Closed System Drug Transfer Device

BD PhaSeal™ Application products	68
BD PhaSeal™ Connector	65, 66
BD PhaSeal™ Injector	64
BD PhaSeal™ Protector	63
PhaSeal™ Administration products	66–67

Closed Systems - Needleless IV Access	
BD Q-Syte™	26
Combined Spinal / Epidural Anaesthesia Sets	57
BD Durasafe Plus!™ with CSE lock	57
Combi plugs	29
Connecta™ 3-way Stopcock	28
Connecta™ Extension Tube	24
Connecta™ Extension Tube for Infusion	24
Connector, BD PhaSeal™	65
Cyto-Ad Z®*	67

D

Discardit™	14
Dressing	43
Durasafe Plus!™	57

E

Eclipse™	8
Emerald™	13
Emerald™ PRO	21
Epidural Products	55, 55–56
Epilor™	56
Extension Tubes	24
E-Z Scrub™	59

F

FloPro™	29
Flushing of Vascular Access Devices	48–49

G

General Conditions of Sales	72
-----------------------------	----

H

Heimlich Chest Drain	59
Hub Cutter	23
Hypodermic Needles	5
BD Microlance™	5

I

Indwelling Cannulae	
BD Angiocath™	42
BD Insyte-N™	41
BD Insyte-N™ Autoguard™	36
BD Insyte-N™ Autoguard™ Winged	36
BD Insyte™ Autoguard™	36
BD Insyte™ Autoguard™ BC	37
BD Insyte™ Autoguard™ Winged	36
BD Insyte-W™	41
BD Neoflon™	43
BD Saf-T-Intima™	38
BD Venflon™	40
BD Venflon™ Pro	39
BD Venflon™ Pro Safety	35
Infusion Accessories	29
Infusion Adapter, BD PhaSeal™	67
Infusion Sets	25
Infusion Therapy	24–31
Injection Systems	5–23
Injector, BD PhaSeal™	64
Insyte-A™	50
Insyte-N™	41
Insyte-N™ Autoguard™	36
Insyte™	41
Insyte™ Autoguard™	36
Insyte™ Autoguard™ BC	37
Insyte™ Autoguard™ with wings	36
Insyte-W™	41
Integra	11
IV Dressing	43

L

Large Size Syringes	17
Loss-of-Resistance Syringe	56
Low Dose Syringes	16

M

Microlance™	5
M/P (Multi Purpose) Procedure Needle	51

1

2

3

4

5

6

7

8

9

10

11

12

13

Index

Alphabetical Index

N

- Needleless IV Access - Closed System ... 26
 - BD Q-Syte™ 26
- Needles and Syringes 13
 - BD Discardit™ 14
 - BD Emerald™ 13
 - BD Large Size Syringes 17
 - BD Plastipak™ 15
 - BD Syringe Luer-Lok™ Tip 15
- Neoflon™ 43
- Nexiva™ Diffusics™ 34
- Nexiva™ Dual Port "Y" with one BD Q-Syte™ 33
- Nexiva™ Dual Port "Y" with two BD Q-Syte™ 33
- Nexiva™ Single Port 33
- Nexiva™ with BD Connecta™ Stopcock .. 33
- Nexiva™ with BD Connecta™ Stopcock and BD Q-Syte™ 33

O

- Obturator 44, 44–45
- One-way Valve 30
- Oral Dispensing Syringes
 - BD Oral Dispensing Syringes 18

P

- Passive Reuse Prevention Syringe 21
 - BD Emerald™ PRO 21
- Percutaneous Injuries 6
- Peripheral Vascular Access Systems 32–47
- Perisafe Plus!™ Epidural Procedure Kit. 55
 - Mini-Perisafe Plus!™ Epidural Procedure Kit 55
 - Perisafe Plus!™ 55
- Persist™ Plus 31
- PhaSeal™
 - BD PhaSeal™ Administration products 66, 67
 - BD PhaSeal™ Application products 68
 - BD PhaSeal™ Connector 65

- BD PhaSeal™ Injector 64
- BD PhaSeal™ Protector 63

- Plastipak™ 15
- Plastipak™ Large Size Syringes 17
- Plugs 29
- PosiFlush™ 48

Q

- Q-Syte™ 26, 27
- Quality Management with BD 74
- Quincke Type Point Spinal Needle... 53–54

R

- Reuse Prevention 19
- Reuse Prevention Syringe
 - BD Emerald™ PRO 21

S

- Safety Box For Disposal of Hypodermic Needles
 - BD™ Hub Cutter 23
- Safety-Engineered Needles and Syringes 8
 - BD Eclipse™ 8
 - BD Integra™ 11
 - BD SafetyGlide™ 9
- SafetyGlide™ 9, 10
- SafetyGlide™ - with BD Luer-Lok™ Syringe packed aside 9
- Saf-T-Intima™ 38
- Secondary Set, BD PhaSeal™ 66
- Sharps Disposal Systems 60–61
- Single Components for Epidural Anaesthesia 56
 - BD Epilor™ 56
 - Tuohy and Tuohy Weiss Epidural Needle 56
- Skin Prep 31
 - BD Persist™ Plus 31
- Social Responsibility - Trusted Partners.. 70
- SoloShot™ Mini 22
- Spinal Anaesthesia 52–54

- Spinal Needle Quincke 53
- Stopcock 28
- Stopcock for IV Therapy and Pressure Monitoring 28
 - BD Connecta™ 3-way Stopcock 28
- Surgical Accessories 58–59
- Symbol Guidance 74
- Syringe Luer-Lok™ Tip 15
- Syringe With Passive Reuse Prevention Technology 20, 21
 - BD Emerald™ PRO 21

T

- Transfer Needle 29
 - BD FloPro™ Transfer Needle 29
- Tuohy and Weiss 56
- Tuohy Epidural Needles 56

U

- Universal Fill Needles 7
 - BD™ Blunt Fill Needle 7

V

- Veca-C™ 43
- Vecafix™ 43
- Venflon™ 40
- Venflon™ Pro 39
- Venflon™ Pro Safety 35
- Vialon 46

W

- Weiss Epidural Needles 56
- Westcott Biopsy Needle 58
- Whitacre Pencil Point Spinal Needle 52

Index

Alphabetical Index

1

2

3

4

5

6

7

8

9

10

11

12

13

Becton Dickinson International –
Branch of Becton, Dickinson B.V.
Erembodegem-Dorp 86
9320 Erembodegem
Belgium
VAT BE 0443 774 010
RPR Dendermonde
Tel: +32 (0)53 720 211
Fax: +32 (0)53 720 200
www.bd.com



Helping all people
live healthy lives

Whilst every effort has been made to ensure the accuracy of the information contained in our product catalogue, the content is naturally subject to change and we cannot guarantee its accuracy. We reserve the right to make changes to the product catalogue without notice. We do not accept any liability for any error or omission contained in our product catalogue.

BD, BD Logo and all other trademarks - unless stated otherwise - are property of Becton, Dickinson and Company. © 2014 BD